

UM-HMC-3B cells | 305715

General information

Description

UM-HMC-3B is a human mucoepidermoid carcinoma cell line established from a primary mucoepidermoid carcinoma of the hard palate salivary glands of a 73-year-old Caucasian female patient. The line is part of the University of Michigan Human Mucoepidermoid Carcinoma (UM-HMC) panel. Like UM-HMC-1, UM-HMC-3B harbors the CRTC1-MAML2 gene fusion (MECT1-MAML2), the oncogenic driver in the majority of mucoepidermoid carcinomas, which results from a t(11;19)(q21;p13) translocation and drives CREB and Notch pathway activation.

UM-HMC-3B is applicable in studies of hard palate and salivary gland mucoepidermoid carcinoma biology, CRTC1-MAML2 signalling, and comparative analyses with other UM-HMC panel lines. Together with UM-HMC-1 and other panel members, UM-HMC-3B enables identification of common and divergent molecular features across mucoepidermoid carcinomas of different anatomical sites, ages, and genders. It is also suitable for drug sensitivity assays targeting fusion oncoproteins and downstream CREB/Notch effectors.

Organism

Human

Tissue

Hard palate, salivary gland

Disease

Hard palate mucoepidermoid carcinoma

Metastatic site

Primary tumor site (hard palate)

Applications

Salivary gland carcinoma research; CRTC1-MAML2 fusion biology; hard palate MEC biology; CREB/Notch pathway; comparative UM-HMC panel studies; drug sensitivity assays

Synonyms

University of Michigan-Human Mucoepidermoid Carcinoma-3B

Characteristics

Age

73 years

Gender

Female

Ethnicity

Caucasian

Morphology

Epithelial-like

Cell type

Epithelial cells

Growth properties

Adherent

UM-HMC-3B cells | 305715

Regulatory Data

Citation	UM-HMC-3B (Cytion catalog number 305715)
Biosafety level	1
NCBI_TaxID	9606
CellosaurusAccession	CVCL_Y472
GMO Status	No genetic modification; somatic CRTC1-MAML2 fusion acquired during tumor development

Biomolecular Data

Mutational profile	Mutation: Gene fusion, CRTC1 + HGNC, MAML2, Name(s)=CRTC1-MAML2, MECT1-MAML2.
--------------------	---

Handling

Culture Medium	DMEM:Ham's F12 (1:1), w: 3.1 g/L Glucose, w: 2.5 mM L-Glutamine, w: 15 mM HEPES, w: 0.5 mM Sodium pyruvate, w: 1.2 g/L NaHCO ₃ (Cytion article number 820400a)
Supplements	Supplement the medium with 10% FBS
Dissociation Reagent	Accutase
Doubling time	approx. 24 to 48 hours
Subculturing	Remove medium, wash with PBS without calcium and magnesium, cover with Accutase, incubate 8–10 min at RT, resuspend in medium, centrifuge 300×g 3 min, discard supernatant, reseed in fresh medium.
Split ratio	1 to 5
Seeding density	2-4*10 ⁴ /cm ²
Fluid renewal	Every 2 to 3 days

UM-HMC-3B cells | 305715

Freeze medium

As a cryopreservation medium, we use complete growth medium (including FBS) + 10% DMSO for adequate post-thaw viability, or CM-1 (Cytion catalog number 800100), which includes optimized osmoprotectants and metabolic stabilizers to enhance recovery and reduce cryo-induced stress.

Thawing and Culturing Cells

1. Confirm that the vial remains deeply frozen upon delivery, as cells are shipped on dry ice to maintain optimal temperatures during transit.
2. Upon receipt, either store the cryovial immediately at temperatures below -150°C to ensure the preservation of cellular integrity, or proceed to step 3 if immediate culturing is required.
3. For immediate culturing, swiftly thaw the vial by immersing it in a 37°C water bath with clean water and an antimicrobial agent, agitating gently for 40-60 seconds until a small ice clump remains.
4. Perform all subsequent steps under sterile conditions in a flow hood, disinfecting the cryovial with 70% ethanol before opening.
5. Carefully open the disinfected vial and transfer the cell suspension into a 15 ml centrifuge tube containing 8 ml of room-temperature culture medium, mixing gently.
6. Centrifuge the mixture at 300 x g for 3 minutes to separate the cells and carefully discard the supernatant containing residual freezing medium.
7. Gently resuspend the cell pellet in 10 ml of fresh culture medium. For adherent cells, divide the suspension between two T25 culture flasks; for suspension cultures, transfer all the medium into one T25 flask to promote effective cell interaction and growth.
8. Adhere to established subculture protocols for continued growth and maintenance of the cell line, ensuring reliable experimental outcomes.

Incubation Atmosphere

37°C, 5% CO₂, humidified atmosphere.

Flask Coating

Use Collagen coated flasks

Shipping Conditions

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately -78 °C throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

UM-HMC-3B cells | 305715

Storage Conditions

For long-term preservation, place vials in vapor-phase liquid nitrogen at about -150 to -196 °C. Storage at -80 °C is acceptable only as a short interim step before transfer to liquid nitrogen.

Quality Control & Molecular Analysis

Sterility

Mycoplasma contamination is excluded using both PCR-based assays and luminescence-based mycoplasma detection methods.

To ensure there is no bacterial, fungal, or yeast contamination, cell cultures are subjected to daily visual inspections.