

NCI-H196 Cells | 300390

Renseignements généraux

Description

NCI-H196 is a small-cell lung cancer (SCLC) cell line used to study the mechanisms of cancer progression, chemotherapy resistance, and cellular responses to oxidative stress. Research involving NCI-H196 has demonstrated its sensitivity to the cytotoxic effects of pyrrolidine dithiocarbamate (PDTC), a pro-oxidant agent. PDTC induces S-phase cell cycle arrest and significantly reduces the viability of NCI-H196 cells in a dose-dependent manner. This cytotoxicity is attributed to the induction of oxidative stress, as evidenced by increased reactive oxygen species (ROS) and changes in the expression of oxidative stress-related genes. The addition of antioxidants like N-acetyl-L-cysteine (NAC) can effectively reverse PDTC-induced cytotoxicity, confirming the role of oxidative stress in cell death.

Further studies have shown that PDTC enhances the cytotoxicity of cisplatin, a first-line chemotherapy drug used for SCLC treatment. Combining low doses of cisplatin with non-toxic concentrations of PDTC leads to synergistic cytotoxicity in NCI-H196 cells. This combination therapy is believed to be effective due to PDTC's downregulation of ATP7A, a copper efflux transporter associated with cisplatin resistance. By inhibiting ATP7A, PDTC may increase intracellular copper and sensitize NCI-H196 cells to cisplatin, highlighting its potential as an adjunct therapy for SCLC.

Organism

Human

Tissue

Lung

Disease

Lung small cell carcinoma

Metastatic site

Pleural effusion

Applications

3D cell culture, Cancer research

Synonyms

NCI-H196, H-196, NCIH196

Caractéristiques

Age

68 years

Gender

Male

Ethnicity

European

Growth properties

Adherent

Données réglementaires

NCI-H196 Cells | 300390

Citation	NCI-H196 (Cytion catalog number 300390)
-----------------	---

Biosafety level	1
------------------------	---

NCBI_TaxID	9606
-------------------	------

CellosaurusAccession	CVCL_1509
-----------------------------	-----------

Données biomoléculaires

Manipulation

Culture Medium	RPMI 1640, w: 2.0 mM stable Glutamine, w: 2.0 g/L NaHCO ₃ (Cytion article number 820700a)
-----------------------	--

Supplements	Supplement the medium with 10% FBS
--------------------	------------------------------------

Dissociation Reagent	Accutase
-----------------------------	----------

Subculturing	Remove the old medium from the adherent cells and wash them with PBS that lacks calcium and magnesium. For T25 flasks, use 3-5 ml of PBS, and for T75 flasks, use 5-10 ml. Then, cover the cells completely with Accutase, using 1-2 ml for T25 flasks and 2.5 ml for T75 flasks. Let the cells incubate at room temperature for 8-10 minutes to detach them. After incubation, gently mix the cells with 10 ml of medium to resuspend them, then centrifuge at 300xg for 3 minutes. Discard the supernatant, resuspend the cells in fresh medium, and transfer them into new flasks that already contain fresh medium.
---------------------	---

Freeze medium	As a cryopreservation medium, we use complete growth medium (including FBS) + 10% DMSO for adequate post-thaw viability, or CM-1 (Cytion catalog number 800100), which includes optimized osmoprotectants and metabolic stabilizers to enhance recovery and reduce cryo-induced stress.
----------------------	---

NCI-H196 Cells | 300390

Thawing and Culturing Cells

1. Confirm that the vial remains deeply frozen upon delivery, as cells are shipped on dry ice to maintain optimal temperatures during transit.
2. Upon receipt, either store the cryovial immediately at temperatures below -150°C to ensure the preservation of cellular integrity, or proceed to step 3 if immediate culturing is required.
3. For immediate culturing, swiftly thaw the vial by immersing it in a 37°C water bath with clean water and an antimicrobial agent, agitating gently for 40-60 seconds until a small ice clump remains.
4. Perform all subsequent steps under sterile conditions in a flow hood, disinfecting the cryovial with 70% ethanol before opening.
5. Carefully open the disinfected vial and transfer the cell suspension into a 15 ml centrifuge tube containing 8 ml of room-temperature culture medium, mixing gently.
6. Centrifuge the mixture at $300 \times g$ for 3 minutes to separate the cells and carefully discard the supernatant containing residual freezing medium.
7. Gently resuspend the cell pellet in 10 ml of fresh culture medium. For adherent cells, divide the suspension between two T25 culture flasks; for suspension cultures, transfer all the medium into one T25 flask to promote effective cell interaction and growth.
8. Adhere to established subculture protocols for continued growth and maintenance of the cell line, ensuring reliable experimental outcomes.

Incubation Atmosphere

37°C , 5% CO_2 , humidified atmosphere.

Shipping Conditions

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately -78°C throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

Storage Conditions

For long-term preservation, place vials in vapor-phase liquid nitrogen at about -150 to -196°C . Storage at -80°C is acceptable only as a short interim step before transfer to liquid nitrogen.

Contrôle de la qualité et analyse moléculaire

NCI-H196 Cells | 300390

Sterility

Mycoplasma contamination is excluded using both PCR-based assays and luminescence-based mycoplasma detection methods.

To ensure there is no bacterial, fungal, or yeast contamination, cell cultures are subjected to daily visual inspections.