

IM-9 Cells | 302151

Renseignements généraux

Description

IM-9 is a human lymphoblastoid cell line established in 1967 from the bone marrow of an adult woman diagnosed with multiple myeloma. Originally believed to be derived from myeloma cells, subsequent research, including findings published by Pellat-Deceunynk et al. in 1995, has indicated that IM-9 cells are more accurately classified as Epstein-Barr Virus-positive (EBV+) B-lymphoblastoid cells rather than malignant myeloma cells. This distinction is crucial for researchers using this cell line, as it affects the interpretation of experimental outcomes related to myeloma studies.

IM-9 cells have been extensively characterized in the literature and are noted for synthesizing Immunoglobulin G (IgG). They are also known to express receptors for insulin and calcitonin, making them valuable for studying hormone-receptor interactions. Additionally, these cells express BCL2 mRNA, a gene involved in regulating apoptosis, which is often studied in the context of cancer and immune cell survival. Due to their high expression of insulin receptors, IM-9 cells are frequently used in research focused on insulin signaling and metabolic disorders, providing insights into insulin resistance mechanisms.

The IM-9 cell line remains a significant resource for various research applications, particularly in the fields of immunology, cancer biology, and metabolic studies. However, given the revised understanding of their origin, it is critical to use IM-9 cells with the recognition that they are not representative of malignant myeloma cells. As always, these cells are intended exclusively for in vitro research and are not suitable for therapeutic or in vivo use.

Organism Human

Tissue Bone marrow

Synonyms IM 9, IM9, GM04680

Caractéristiques

Age Unspecified

Gender Female

Ethnicity Caucasian

Morphology Round cells in cluster

Cell type B lymphoblast

Growth properties Suspension

IM-9 Cells | 302151

Données réglementaires

Citation	IM-9 (Cytion catalog number 302151)
Biosafety level	2
NCBI_TaxID	9606
CellosaurusAccession	CVCL_1305

Données biomoléculaires

Antigen expression	CD19+, CD20+, CD23+, CD27+, CD80+, CD83+, CD138+, MHC I+, MHC II+
Viruses	EBV+ free of human pathogenic viruses SV40, JC/BK, HBV, HCV, HIV.

Manipulation

Culture Medium	RPMI 1640, w: 2.0 mM stable Glutamine, w: 2.0 g/L NaHCO ₃ (Cytion article number 820700a)
Supplements	Supplement the medium with 10% heat-inactivated FBS
Subculturing	Gently homogenize the cell suspension in the flask by pipetting up and down, then take a representative sample to determine the cell density per ml. Dilute the suspension to achieve a cell concentration of 1×10^5 cells/ml with fresh culture medium, and aliquot the adjusted suspension into new flasks for further cultivation.
Fluid renewal	2 to 3 times per week
Post-Thaw Recovery	Fast
Freeze medium	As a cryopreservation medium, we use complete growth medium (including FBS) + 10% DMSO for adequate post-thaw viability, or CM-1 (Cytion catalog number 800100), which includes optimized osmoprotectants and metabolic stabilizers to enhance recovery and reduce cryo-induced stress.

IM-9 Cells | 302151

Thawing and Culturing Cells

1. Confirm that the vial remains deeply frozen upon delivery, as cells are shipped on dry ice to maintain optimal temperatures during transit.
2. Upon receipt, either store the cryovial immediately at temperatures below -150°C to ensure the preservation of cellular integrity, or proceed to step 3 if immediate culturing is required.
3. For immediate culturing, swiftly thaw the vial by immersing it in a 37°C water bath with clean water and an antimicrobial agent, agitating gently for 40-60 seconds until a small ice clump remains.
4. Perform all subsequent steps under sterile conditions in a flow hood, disinfecting the cryovial with 70% ethanol before opening.
5. Carefully open the disinfected vial and transfer the cell suspension into a 15 ml centrifuge tube containing 8 ml of room-temperature culture medium, mixing gently.
6. Centrifuge the mixture at $300 \times g$ for 3 minutes to separate the cells and carefully discard the supernatant containing residual freezing medium.
7. Gently resuspend the cell pellet in 10 ml of fresh culture medium. For adherent cells, divide the suspension between two T25 culture flasks; for suspension cultures, transfer all the medium into one T25 flask to promote effective cell interaction and growth.
8. Adhere to established subculture protocols for continued growth and maintenance of the cell line, ensuring reliable experimental outcomes.

Incubation Atmosphere

37°C , 5% CO_2 , humidified atmosphere.

Shipping Conditions

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately -78°C throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

Storage Conditions

For long-term preservation, place vials in vapor-phase liquid nitrogen at about -150 to -196°C . Storage at -80°C is acceptable only as a short interim step before transfer to liquid nitrogen.

Contrôle de la qualité et analyse moléculaire

IM-9 Cells | 302151

Sterility

Mycoplasma contamination is excluded using both PCR-based assays and luminescence-based mycoplasma detection methods.

To ensure there is no bacterial, fungal, or yeast contamination, cell cultures are subjected to daily visual inspections.