

[illegible]

**Age**                      00000 000000

**Citation** BJ (XXXX XXXXXX) Cytion 305222)

BJ  | 305222

  - 

**Karyotype**

BJ                       



**Culture Medium**

DMEM, w: 4.5 g/L  w: 4 mM L- w: 3.7 g/L NaHCO<sub>3</sub>, w: 1.0 mM  (  Cytion 820300a)

**Supplements**

 10% FBS, 20 ng/mL bFGF

**Dissociation Reagent**



**Subculturing**

**Freeze medium**

**Thawing and Culturing Cells**

1.                       
2.                       
3.                       
4.                       
5.                       
6.                       
7.                       
8.                       

**Incubation Atmosphere**

37°C, 5% CO<sub>2</sub>, 

BJ  | 305222

Flask Coating 


Freezing  
Procedure

  -78°C

Shipping  
Conditions

  -78°C

Storage  
Conditions

   -150  196  

Sterility

    PCR  