

XXXXXXXXXXXXXXXXXXXX

XXXXXXXXXX XXXXXXXXXXXXXXXXXXXXXXXX

Citation	RWPE-1 (GenBank: F05217)
Biosafety level	RWPE-1 (GenBank: F05217) is a BSL-1/2 (BSL-1/2) cell line, derived from a human embryo. It is a BSL-1/2 cell line, derived from a human embryo. It is a BSL-1/2 cell line, derived from a human embryo.
NCBI_TaxID	9606
CellosaurusAccession	CVCL_3791

~~XXXXXXXXXX~~ ~~XXXX~~- ~~XXXXXXXXXXXXXXXXXXXX~~

XXXXXX

Supplements [XXXXX XXXXXXXX 0.05 μg/μl BPE, 5 μg/μl EGF. XXX XXXXX XX XXXXXXX XXXXXXX. XXXX BPE X-EGF X-10 μg/μl, XXXXXXX XXXXXXX XXXXXXX, XXX](#)

Subculturing

Thawing and Culturing Cells

1. 取出冻存管，将冻存管放入37℃水浴中快速解冻，取出后立即放入冰浴中冷却，离心（1000g，5min），弃去上清液，加入1ml新鲜培养基，重悬细胞。
2. 将细胞悬液接种到含有10ml新鲜培养基的25cm²培养瓶中，置于37℃、5% CO₂培养箱中培养。
3. 待细胞贴壁后，将培养基更换为新鲜培养基，继续培养。
4. 待细胞长满瓶底后，将细胞悬液接种到含有10ml新鲜培养基的25cm²培养瓶中，置于37℃、5% CO₂培养箱中培养。
5. 待细胞长满瓶底后，将细胞悬液接种到含有10ml新鲜培养基的25cm²培养瓶中，置于37℃、5% CO₂培养箱中培养。
6. 待细胞长满瓶底后，将细胞悬液接种到含有10ml新鲜培养基的25cm²培养瓶中，置于37℃、5% CO₂培养箱中培养。
7. 待细胞长满瓶底后，将细胞悬液接种到含有10ml新鲜培养基的25cm²培养瓶中，置于37℃、5% CO₂培养箱中培养。
8. 待细胞长满瓶底后，将细胞悬液接种到含有10ml新鲜培养基的25cm²培养瓶中，置于37℃、5% CO₂培养箱中培养。

Incubation Atmosphere	37°C, 5% CO ₂ ,                                   

