

MM.1S | 305304

Cell Line

Products IgA

Mutational profile KRAS, p.Gly12Ala (c.35G>C), TRAF3, p.Val536_Asn545delValPheValAlaGlnThrValLeuGluAsninsA 1630delTCTTTGTGGCCAAACTGTTCTAGAAA,

Media

Culture Medium RPMI 1640, w: 2.0 mM, w: 2.0 g/L NaHCO3 (Cytion 820700a)

Supplements 10% FBS

Dissociation Reagent

Subculturing 15 min PBS

Freeze medium (FBS) + 10% DMSO

Thawing and Culturing Cells

- 1. ...
- 2. ... -150°C ...
- 3. ... 37 ...
- 4. ... 70% ...
- 5. ... 15 min 8 min ...
- 6. ... -300 x g 3 ...
- 7. ... -10 min ...
- 8. ...

