

XXXXXXXXXXXX

XXXXXXXXXX XIXX-XXXXXXXXXXXXXXXXXXXX

GES-1 | 305428

Shipping
Conditions

GES-1 is a dry powder reagent. It is stable at room temperature for up to 12 months. It is shipped in a dry ice container at -78°C.

Storage
Conditions

GES-1 should be stored at -150 to -196 °C. It is stable for up to 12 months.

GES-1 is a dry powder reagent. It is stable at room temperature for up to 12 months.

Sterility

GES-1 is sterile and ready to use. It is not recommended to use GES-1 for PCR. It is recommended to use GES-1 for sequencing.