

DOHH-2 | 305537

DOHH-2

Description	DOHH-2 is a human B cell line (DLBCL) derived from a patient with diffuse large B-cell lymphoma (DLBCL). It is a cell line that is used for the study of B cell biology and for the development of new therapies. DOHH-2 is a cell line that is used for the study of B cell biology and for the development of new therapies. DOHH-2 is a cell line that is used for the study of B cell biology and for the development of new therapies.
-------------	---

Organism	Human
----------	-------

Tissue	Diffuse large B-cell lymphoma
--------	-------------------------------

Disease	Diffuse large B-cell lymphoma
---------	-------------------------------

Synonyms	DOHH2, DoHH-2, DOHH-2
----------	-----------------------

DOHH-2

Age	60 years
-----	----------

Gender	Male
--------	------

Ethnicity	Chinese
-----------	---------

Growth properties	Adherent
-------------------	----------

DOHH-2

Citation	DOHH-2 (Cytion 305537)
----------	------------------------

Biosafety level	1
-----------------	---

NCBI_TaxID	9606
------------	------

CellosaurusAccession	CVCL_1179
----------------------	-----------

GMO Status	GMO-S1: DOHH-2 (Cytion 305537) is a cell line derived from a patient with diffuse large B-cell lymphoma (DLBCL). It is a cell line that is used for the study of B cell biology and for the development of new therapies.
------------	---

DOHH-2

DOHH-2 | 305537

Viruses

Epstein-Barr virus (EBV).

Cell line

Culture Medium

RPMI 1640, w: 2.0 mM NaCl , w: 2.0 g/L NaHCO_3 (Cytion 820700a)

Supplements

10% FBS

Subculturing

Cells are grown in 15 ml flasks in the presence of 10% FBS. Cells are harvested by trypsinization (3-5 min) and resuspended in PBS.

Seeding density

$3-5 \times 10^5/\text{ml}$

Freeze medium

Cells are frozen in the presence of 10% FBS + 10% DMSO.

Thawing and Culturing Cells

1. Cells are thawed in a water bath at 37°C.
2. Cells are washed with PBS and resuspended in 10% FBS.
3. Cells are seeded in 15 ml flasks.
4. Cells are grown in the presence of 10% FBS.
5. Cells are harvested by trypsinization (3-5 min) and resuspended in PBS.
6. Cells are counted and seeded in 15 ml flasks.
7. Cells are grown in the presence of 10% FBS.
8. Cells are harvested by trypsinization (3-5 min) and resuspended in PBS.

Incubation Atmosphere

37°C, 5% CO_2

Shipping Conditions

Cells are shipped in the presence of 10% FBS + 10% DMSO.

DOHH-2 | 305537

Storage Conditions

XXXXXXXX XXXX XXX, XX XXXXX XX XXXXXXXXXXX XXXXX XXXXX XXXX XX XXX XXXXXXXXXXX XX X-150 XX 196 XXXXX XXXXXXX. XXXXX XXXXXXX

Sterility

XXXXXXXXXXXXXXXXXXXX XXXX XXXXXXXX XXXXXXXX PCR XXXXXXXX XXXXXXXX XXXXXXXXXXXXXXXXXXXX

[illegible]