

XXXX XXXXXX | 305351

XXXXX XXXXXX

Description	XX XXXXX "XXXXX" XX XX XXXX XX XXXXXXX XXXXXXX XX XX B XXXXXXX (DLBCL), XXXXX XXXXX XX XXXXX XX XXXX X-DLBCL, XX XX XX XXXXXXX XXXXXXX XX XX XXXXX XXXXX XX XXXXX XXXXX XX XXXXX XXXXXXXXX XX XXXXXXXXXX (Ig), XX XXXXX XXXX XX
Organism	XX
Tissue	XX XXXXX
Disease	XXXXXXXX XXXXXXX XX XX B XXXXXXX XXXX XX B XXXX XXXXXXX
Applications	XXXXX XXXX XX-XXXX, XXXXXXXXXX
Synonyms	XXXXX

XXXXXXXXXX

Age	XXXXX
Gender	XXXX
Ethnicity	XXXXXXXX
Morphology	XXXXXXXXXX
Growth properties	XXXXX

XXXXXXXXX XXXXXXXXXXXXXXXX

Citation	XXXXX (XXXX XXXXX Cytion 305351)
Biosafety level	1
NCBI_TaxID	9606
CellosaurusAccession	CVCL_3611

XXXX XXXXXX | 305351

XIX-XI-XXXXXXXXXXXX

Antigen expression	CD10+, CD19+, CD20+, CD38+, CD23-, CD39-
Viruses	EBV-, HBV-, HCV-, HIV-1-, HIV-2-, HTLV-1/2-, MLV-, SMRV-
Mutational profile	XXXXXXXX KRAS, XXXXXX p.Gly13Asp (c.38G>A), XXXXXXXXXXXX, XXXXXXXX TP53, XXXXXX p.Glu198fs (c.592_596del5), XX XXXXX

XXXXXX

Culture Medium	RPMI 1640, w: 2.0 mM CaCl_2 , w: 2.0 g/L NaHCO_3 (5% CaCl_2 5% Cytion 820700a)
Supplements	10% FBS
Dissociation Reagent	Trypsin
Doubling time	50 days
Subculturing	5-6 days, 1:5 x 10
Seeding density	2×10^4 cells/cm ²
Fluid renewal	2-3 days
Freeze medium	10% FBS + 10% DMSO

XXXX XXXXXX | 305351

**Thawing and
Culturing Cells**

1. XXX XXXXXXXXXX XXXX XXXX XXXXXXXXXX XXX XXXXXXX, XXX XXXXXXX XXXXXXX XXXX XXX XXXXXXX XXX XXXXXXXXXX XXXXXXXXXX XXXXXXXXXX
2. XXX XXXX XXXXXXX, XXXX XXX XXXXXXXXXX XXXXXXXXXX XXX XXXXXXXXXX XXXXXXX XXX -150°C XXX XXXXXXX XXX XXXXXXX XXXXXXX XXXXXXX, XXX XXXXXXX
3. XXXXXXX XXXXXXX XXXXXXX, XXXXXXX XXX XXXXXXXXXX XXXXXXXXXX XXX XXX XXXXXXX XXXXXXX XXX XXXXXXXXXX XXX 37 XXXXXXX XXXXXXXXXX XXX XXX XXXXXXX
4. XXX XXX XXX XXXXXXX XXXXXXX XXXXXXXXXX XXXXXXX, XXXX XXX XXXXXXXXXX XXXXXXXXXX XXXXXXX 70% XXXX XXXXXXXXXX
5. XXX XXXXXXXXXX XXX XXXXXXXXXX XXXXXXXXXX XXXXXXX XXX XXXXXXX XXXXXXX XXXXXXXXXX XXXXXXXXXX XXXX 15 μ" XXXX 8 μ" XXX XXX XXXXXXX XXXXXXX
6. XXXXXXX XXX XXXXXXXXXX XXXXXXXXXX XXX-300 x g XXXX 3 XXXX XXX XXXXXXXXXX XXX XXXXXXX, XXXXXXX XXXXXXXXXX XXX XXXXXXX XXXXXXXXXX XXXXXXX XXX
7. XXXXXXX XXXXXXXXXX XXX XXXX XXXXXXX XXX-10 μ" XXX XXX XXXXXXX XXXX. XXXX XXXX XXXXXXX, XXXX XXX XXXXXXXXXX XXX XXX XXXXXXXXXX XXXXXXX T2
8. XXXXXXXXXX XXX XXXXXXXXXX XXX-XXXXXXX XXXXXXXXXX XXXXXXX XXXXXXXXXX XXX XXX XXXXXXX, XXX XXXXXXXXXX XXXXXXXXXX XXXXXXXXXX XXXXXXXXXX

**Incubation
Atmosphere**

37°C, 5%_{CO2}, XXXXXXXXXX XXXX

**Shipping
Conditions**

XXXXXXXX XXXX XXXXXXXXXX XXXXXXXXXX XXXX XXX XXXXXXXXXX XXXXXXXXXX XXXXXXXXXX XXX XXXX XXXXXXX XXX XXXXXXX XXX XXXXXXX XXX XXXXXXXXXX XXX -78°C

**Storage
Conditions**

XXXXXXXX XXXXXXX XXXX, XXX XXXXXXX XXX XXXXXXXXXX XXXXXXX XXXXXXX XXXX XXX XXX XXXXXXXXXX XXX -150 XXX 196 XXXXXXX XXXXXXXXXX. XXXXXXX XXXXXXX

XXXXX XXXXXXX XXXXXXXXXX XXXXXXXXXX

Sterility

XXXXXXXX XXXXXXXXXX XXXX XXXXXXXXXX XXXXXXXXXX XXXXXXXXXX PCR XXXXXXXXXX XXXXXXX XXXXXXXXXX XXXXXXXXXX XXXXXXXXXX XXXXXXXXXX

XXX XXXXXXX XXXX XXXXXXX XXXXXXX, XXXXXXXXXX XXX XXXX, XXXXXXXXXX XXXXXXX XXXXXXX XXXXXXXXXX XXXXXXXXXX XXXXXXXXXX