

Product sheet

U-CH1 | 305885

Cell line: U-CH1

Mutational profile TP53, p.Pro72Arg (c.215C>G),

Cell line: U-CH1

Culture Medium IMDM, w: 4.5 g/L, w: 4 mM L- , w: 25 mM HEPES, w: 1.0 mM , w: 3.024 g/L NaHCO3 (Cytion)

Supplements 10% FBS

Dissociation Reagent

Doubling time ~

Fluid renewal 2 x 3

Freeze medium (FBS) + 10% DMSO

- Thawing and Culturing Cells**
- 1.
 - 2.
 - 3.
 - 4.
 - 5.
 - 6.
 - 7.
 - 8.

U-CH1 | 305885

Incubation Atmosphere 37°C, 5% CO₂, humidified

Flask Coating Cell binding

Shipping Conditions The cells are shipped in a dry ice container at -78°C. The cells are shipped in a dry ice container at -78°C.

Storage Conditions The cells are stored at -150 to -196°C in liquid nitrogen. The cells are stored at -150 to -196°C in liquid nitrogen.

Cell line: U-CH1, Cell number: 305885

Sterility The cells are sterile and free of mycoplasmas. The cells are sterile and free of mycoplasmas. The cells are sterile and free of mycoplasmas. The cells are sterile and free of mycoplasmas. The cells are sterile and free of mycoplasmas.