

Product sheet

SU-DHL-8 | 305877

Cell line **Cell line**

Antigen expression Ig+; IgM-, IgG-, IgA-, IgD-, Lambda-, Kappa-

Mutational profile **Cell line** KMT2D, **Cell line** p.Pro648Thrfs*2 (c.1940dupC) (c.1940_1941insC), **Cell line** (Cosmic-CLP=1331038), TP53, **Cell line** p.T (c.700T>A), **Cell line** (Cosmic-CLP=1331038), TP53, **Cell line** p.Arg249Gly (c.745A>G), **Cell line** (Cosmic-CLP=1331038)

Cell line

Culture Medium RPMI 1640, w: 2.0 mM **Cell line**, w: 2.0 g/L NaHCO3 (**Cell line** Cytion 820700a)

Supplements **Cell line** 10% FBS

Dissociation Reagent **Cell line**

Doubling time ~48-72 **Cell line**

Seeding density 0.3-0.5 x 10⁶ **Cell line**

Fluid renewal 2 x 3 **Cell line**

Freeze medium **Cell line**, **Cell line** (FBS) + 10% DMSO **Cell line**

SU-DHL-8 | 305877

**Thawing and
Culturing Cells**

1. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，避免产生气泡。将细胞培养瓶置于37°C孵育箱中，静置15分钟。
2. 将细胞培养瓶置于37°C孵育箱中，静置15分钟。将细胞培养瓶置于37°C孵育箱中，静置15分钟。
3. 将细胞培养瓶置于37°C孵育箱中，静置15分钟。将细胞培养瓶置于37°C孵育箱中，静置15分钟。
4. 将细胞培养瓶置于37°C孵育箱中，静置15分钟。将细胞培养瓶置于37°C孵育箱中，静置15分钟。
5. 将细胞培养瓶置于37°C孵育箱中，静置15分钟。将细胞培养瓶置于37°C孵育箱中，静置15分钟。
6. 将细胞培养瓶置于37°C孵育箱中，静置15分钟。将细胞培养瓶置于37°C孵育箱中，静置15分钟。
7. 将细胞培养瓶置于37°C孵育箱中，静置15分钟。将细胞培养瓶置于37°C孵育箱中，静置15分钟。
8. 将细胞培养瓶置于37°C孵育箱中，静置15分钟。将细胞培养瓶置于37°C孵育箱中，静置15分钟。

**Incubation
Atmosphere**

37°C, 5% CO₂, 饱和湿度

Flask Coating

无涂层

**Shipping
Conditions**

细胞悬液应储存在-150°C至-196°C的液氮中，运输过程中应保持在-78°C以下。

**Storage
Conditions**

细胞悬液应储存在-150°C至-196°C的液氮中，运输过程中应保持在-78°C以下。

细胞悬液 **细胞悬液** **细胞悬液** **细胞悬液**

Sterility

细胞悬液应储存在-150°C至-196°C的液氮中，运输过程中应保持在-78°C以下。