

XXXXXXXXXX XXXXXXXX XXXXXXXX | 300692

XXXXX XXXXX

Description	XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXXX (HEKs) XX XXXX XXXXXXX XXXXXXXXXXX XXXXXXXXXXX XXXXXXXXXXX XXXX XXXX, XXXXXXX XXXX XXX XXXXXXX XXXX HEK XXXXXXX XX XXXXXXXXXXX XXXXXXXXXXX XXXX XX XXXXXXXXXXX XXXXXXX, XX XXXXXXX XXXXXXX XXXX XXXXXXXXXXX XXXXXXXXXXX
-------------	--

Organism	XXXX
----------	------

Tissue	XXXX, XXXXXXXX
--------	----------------

Disease	XXXXX
---------	-------

Applications	XXXXXXXXXXXX, XXXXXXX XXXXXXX, XXXX XXXX, XXXXXXX XXXXXXX UV, XXXXXXXXXXX, XXXXXXX, XXXXXXX XXXXXXX, XXXXXXX XXXXXXX XXXX, XXXXXXX XXXX, XXXXXXX XXXX
--------------	---

XXXXXXXXXXXX

Age	XXXXXXXX
-----	----------

Gender	XXXXXXXX XXXXXXX
--------	------------------

Ethnicity	XXXXXXXX XXXXXXX
-----------	------------------

Morphology	XXXXX XX XXXX XXXXXXX, XXXXXXX XXXXXXXX, XX XXXXXXXX, XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX, XXXXXXX XX 80% XXXXXXX XXXX XXXXXXX XX XXXX
------------	--

Cell type	XXXXXXXXXXXX
-----------	--------------

Growth properties	XXXX
-------------------	------

XXXXXXXXXX XXXXXXXXXXXXXXXX

Citation	XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXXX (XXXX XXXXXXX Cytion 300692)
----------	---

Biosafety level	1
-----------------	---

NCBI_TaxID	9606
------------	------

XXXXXXXXXX XXXX-XXXXXXXXXXXXX

XXXXXXXXXX XXXXXXXXXXXX XXXXXXXX | 300692

XXXXXXXXX

Freeze medium

XXXXXXXXX XXXXXXXX XXXXXXXX, XXXX XXXXXXXX XXXXXXXX XXXXXXXX XXXX (XXXXX FBS) + 10% DMSO XXXX XXXXXXXX XXXXXXXX XXXXXXXX XXXXXXXX XXXXXXXX, XXXX C

Thawing and Culturing Cells

1. XXXX XXXXXXXX XXXX XXXX XXXXXXXX XXXX XXXXXXXX, XXXX XXXXXXXX XXXXXXXX XXXX XXXX XXXX XXXX XXXX XXXXXXXX XXXXXXXX XXXXXXXX XXXXXXXX
2. XXX XXXX XXXXXXXX, XXXX XXX XXXXXXXX XXXXXXXX XXXX XXXXXXXX XXXXXXXX -150°C XXX XXXXXXXX XXX XXXXXXX XXXXXXX XXXXXXX, XXX XXXX
3. XXXXXXX XXXXXXX XXXXXXX, XXXXXXX XXX XXXXXXXX XXXXXXXX XXX XXX XXXXXXX XXXXXXX XXX XXXXXXXX XXX 37 XXXXXXX XXXXXXX XXX XXX XXXX
4. XXXX XXX XXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX, XXXX XXX XXXXXXX XXXXXXX XXXXXXX 70% XXXX XXXXXXX.
5. XXXX XXXXXXX XXX XXXXXXXX XXXXXXX XXXXXXX XXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXX 15 μl XXXX 8 μl XXX XXXX XXXXXXX XXXX
6. XXXXXXX XXX XXXXXXXX XXXXXXX XXXXXXX -300 x g XXXX 3 XXXX XXX XXXXXXX XXX XXXXXXX, XXXXXXX XXXXXXX XXX XXXXXXX XXXXXXX XXXXXXX XXXX
7. XXXXXXX XXXXXXX XXX XXXX XXXXXXX -10 μl XXX XXX XXXX XXXX. XXXX XXXX XXXXXXX, XXXX XXX XXXXXXX XXXX XXX XXXXXXX XXXXXXX T2
8. XXXXXXX XXX XXXXXXXX XXX-XXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXX XXXXXXX, XXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX.

Incubation Atmosphere

37°C, 5% CO₂, XXXXXXX XXXX

Flask Coating

XXX XXXX

Shipping Conditions

XXXXXXXXX XXXX XXXXXXXX XXXXXXX XXXX XXX XXXXXXX XXXXXXX XXXXXXX XXX XXXX XXXXXXX XXX XXXXXXX XXX XXXXXXX XXX XXXXXXX XXX XXXXXXX XXX XXXXXXX XXX -78°C

Storage Conditions

XXXXXXXXX XXXXXXX XXXX, XXX XXXXXXX XXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXX XXX XXX XXXXXXX XXX -150 XXX 196 XXXXXXX XXXXXXX XXXXXXX XXXXXXX

XXXXXXXXX XXXXXXXX XXXXXXXX XXXXXXXX XXXXXXXX

Sterility

XXXXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX PCR XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX

XXXX XXXXXXX XXXX XXXXXXX XXXXXXX, XXXXXXX XXX XXXX, XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX