

CHO-PD-L1 | 305975

CHO-PD-L1

Description	CHO-PD-L1 is a CHO cell line expressing PD-L1 (CD274/B7-H1). The cells are grown in DMEM supplemented with 10% FBS. The cells are used for studying the expression and function of PD-L1 in cancer cells.
-------------	--

Organism	CHO
----------	-----

Tissue	CHO
--------	-----

Disease	PD-L1 (CD274/B7-H1)
---------	---------------------

Applications	CHO-PD-L1 cells are used for studying the expression and function of PD-L1 in cancer cells.
--------------	---

CHO-PD-L1

Age	CHO
-----	-----

Gender	CHO
--------	-----

Morphology	CHO
------------	-----

Cell type	CHO
-----------	-----

Growth properties	CHO
-------------------	-----

CHO-PD-L1

Citation	CHO-PD-L1 (CHO cell line) Cytion: 305975
----------	--

Biosafety level	1
-----------------	---

NCBI_TaxID	10029
------------	-------

CHO-PD-L1 | 305975

**Thawing and
Culturing Cells**

1. 将冻存管放入37°C水浴中快速解冻，解冻后将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。
2. 将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。
3. 将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。
4. 将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。
5. 将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。
6. 将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。
7. 将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。
8. 将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。将细胞悬液转移到含有10 mL完全培养基的细胞培养瓶中，轻轻混匀。

**Incubation
Atmosphere**

37°C, 5% CO₂, 饱和湿度

**Shipping
Conditions**

细胞在含有完全培养基的细胞培养瓶中，在4°C下运输，保质期为7天。

**Storage
Conditions**

细胞在含有完全培养基的细胞培养瓶中，在-150°C下冷冻，保质期为196天。

CHO-PD-L1 | 305975

Sterility

细胞在含有完全培养基的细胞培养瓶中，在4°C下运输，保质期为7天。