

HEK293-Rpn11-HTBH | 305719

HEK293-Rpn11-HTBH

Description	HEK293 - Rpn11-HTBH is a HEK293 (Human Embryonic Kidney 293), JAMM Zn ²⁺ , 26S, (UPS), Rpn11/PSMD14
-------------	--

Organism	
----------	--

Tissue	
--------	--

Disease	HEK293; Rpn11-HTBH
---------	--------------------

Applications	26S; Rpn11/PSMD14; UPS; TEV
--------------	-----------------------------

HEK293-Rpn11-HTBH

Morphology	
------------	--

Cell type	
-----------	--

Growth properties	
-------------------	--

HEK293-Rpn11-HTBH

Citation	HEK293 - Rpn11-HTBH (Cytion 305719)
----------	-------------------------------------

Biosafety level	1
-----------------	---

NCBI_TaxID	9606
------------	------

GMO Status	GMO-S1: HEK293 Rpn11-HTBH (Rpn11/PSMD14 TEV)
------------	--

HEK293-Rpn11-HTBH

HEK293

Culture Medium	DMEM, w: 4.5 g/L, w: 4 mM L-, w: 3.7 g/L NaHCO ₃ , w: 1.0 mM (Cytion 820300a)
----------------	--

Thawing and Culturing Cells

1. 取出冻存管，将冻存管放入37℃水浴中快速解冻，取出后立即放入冰浴中冷却，离心（1000g，5min），弃去上清液，加入1ml新鲜培养基，重悬细胞，接种到24孔板中。
2. 将细胞放入37℃培养箱中培养，待细胞贴壁后，更换新鲜培养基。若细胞未贴壁，可在培养箱中静置过夜，次日再换培养基。
3. 待细胞贴壁后，加入10% FBS培养基，继续培养24-48小时，使细胞进入对数生长期。
4. 待细胞进入对数生长期后，将细胞接种到96孔板中，进行MTT实验。
5. 将细胞接种到96孔板中，每孔接种15-20 μl，接种后在37℃培养箱中培养24-48小时。
6. 待细胞进入对数生长期后，将细胞接种到96孔板中，进行MTT实验。
7. 将细胞接种到96孔板中，进行MTT实验。

Shipping Conditions                  

HEK293-Rpn11-HTBH | 305719

Storage Conditions

XXXXXXXX XXXX XXXX, XX XXXXX XX XXXXXXXXXXX XXXXX XXXXX XXXX XX XXX XXXXXXXXXXX XXX X-150 XX 196 XXXXX XXXXXXX. XXXXX XXXXXXX

