

XXXXXX XXXXXX

XXXXXXXXXX

XIX XIX XIX XIX XIX XIX XIX XIX

Product sheet

HEK293T C6-Luc | 305620

HEK293T C6-Luc - HEK293T cells

Antigen expression	Luc2 (intracellular, extracellular)
Tumorigenic	Yes, HEK293T cells are tumorigenic in SCID-hu mice. HEK293T cells are derived from HEK293 cells, which are derived from HEK293 cells.
Products	HEK293T cells (Luc2, pGL4.50)

HEK293T

Culture Medium	RPMI 1640, w: 2.1 mM NaH_2PO_4 , w: 2.0 g/L NaHCO_3 (HEK293T Cytion 820700a)
Supplements	10% FBS
Dissociation Reagent	Trypsin
Doubling time	24-48 hours
Subculturing	HEK293T cells are typically subcultured at a density of 1×10^4 to 1×10^5 cells per well in 96-well plates. The cells are typically subcultured at a density of 1×10^4 to 1×10^5 cells per well in 96-well plates. The cells are typically subcultured at a density of 1×10^4 to 1×10^5 cells per well in 96-well plates.
Seeding density	$1-3 \times 10^4$ cells/well
Freeze medium	HEK293T cells are typically frozen in a medium containing 10% DMSO and 90% FBS.

