

XXXXXX XXXXXX

[illegible][illegible]ARPE19,      -19, NTC-200, NTC200

XXXXXXXXXXXX

19 □□□□



XXXXXX



ARPE-19 (Cytion 305025)

1

9606

CVCL_0145

XXXXXXXXXX XXXX-XXXXXXXXXXXX

□□□□ □□□□□□ □-Rpe Cralbp □-Rpe-65

ARPE-19 | 305025

Incubation Atmosphere 37°C, 5% CO₂, humidified

Flask Coating none

Freezing Procedure Cells are seeded into 25 cm² flasks and grown to confluence. Cells are then trypsinized and resuspended in freezing medium. The suspension is then aliquoted into 1 mL vials and frozen at -78°C.

Shipping Conditions Cells are shipped on dry ice at -78°C.

Storage Conditions Cells are stored at -150 to -196 °C in liquid nitrogen.

ARPE-19 is a human embryonic kidney cell line.

Sterility

ARPE-19 cells are tested for mycoplasma contamination using PCR. The results of the PCR test are as follows:

ARPE-19 cells are free of mycoplasma contamination.