

FS-C57BL | 400420

General information

|             |                                                                                                                                                                                                                                                                                                                |
|-------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | FS-C57BL is a mouse strain derived from C57BL, which is a highly inbred strain of mouse. FS-C57BL is a highly inbred strain of mouse, derived from C57BL, which is a highly inbred strain of mouse. FS-C57BL is a highly inbred strain of mouse, derived from C57BL, which is a highly inbred strain of mouse. |
| Organism    | Mouse                                                                                                                                                                                                                                                                                                          |
| Tissue      | Whole body                                                                                                                                                                                                                                                                                                     |
| Disease     | None                                                                                                                                                                                                                                                                                                           |

Genetic information

|                   |            |
|-------------------|------------|
| Breed/Subspecies  | C57BL/6J   |
| Gender            | Male       |
| Cell type         | Whole body |
| Growth properties | None       |

Identification and documentation

|                      |                          |
|----------------------|--------------------------|
| Citation             | FS-C57BL (Cytion 400420) |
| Biosafety level      | 1                        |
| NCBI_TaxID           | 10090                    |
| CellosaurusAccession | CVCL_5756                |

Genetic information and documentation

Genetic information

|                |                                                        |
|----------------|--------------------------------------------------------|
| Culture Medium | RPMI 1640, w: 2.0 mM $\text{NaHCO}_3$ (Cytion 820700a) |
|----------------|--------------------------------------------------------|



FS-C57BL | 400420

Freezing  
Procedure

Freeze the sample in a cryovial in a controlled freezing environment at -78°C.

Shipping  
Conditions

Store the sample at -78°C during shipping.

Storage  
Conditions

Store the sample at -150 to -196°C in a cryovial.

FS-C57BL | 400420

Sterility

The sample is sterile and suitable for PCR analysis. The sample is provided in a sterile cryovial.