



XXXXXXXXXXXXXXXXXXXX

**XIX XIX XIX XIX XIX XIX XIX XIX XIX XIX**

**CellosaurusAccession** CVCL E534

XXXXXXXXXX XIXX-XXXXXXXXXXXXXXXXXXXX

XXXXXX XX XXXX EBV

**XXXXXXXXXX**

RPMI 1640, w: 2.0 mM  $\beta$ -mercaptoethanol, w: 2.0 g/L NaHCO<sub>3</sub> (Sigma), 10% Cytion 820700a)

10% FBS 10% FBS

| Age Group | Percentage |
|-----------|------------|
| 18-24     | 15%        |
| 25-34     | 15%        |
| 35-44     | 15%        |
| 45-54     | 15%        |
| 55-64     | 15%        |
| 65-74     | 15%        |
| 75-84     | 15%        |
| 85+       | 10%        |

$$5 \times 10^5 \frac{\text{m}^2}{\text{s}}$$

5 7 10 24

[illegible][illegible]

COX | 302138

**Incubation Atmosphere** 37°C, 5%  $\text{CO}_2$ , humidified

**Flask Coating** Cell binding

**Freezing Procedure** Cells are harvested by centrifugation at 400g for 5 min. The supernatant is removed and the cell pellet is resuspended in freezing medium. The cells are then frozen in a controlled rate freezer and stored at -78°C.

**Shipping Conditions** Cells are shipped in dry ice at -78°C.

**Storage Conditions** Cells are stored at -150 to -196 °C in liquid nitrogen. The storage time is unlimited.

COX | 302138 | 302138

**Sterility** Cells are supplied as a suspension in DMEM/F12 medium. The medium is supplemented with 10% FBS. The cells are tested for mycoplasma contamination by PCR. The results are negative.

Cells are supplied as a suspension in DMEM/F12 medium, supplemented with 10% FBS. The cells are tested for mycoplasma contamination by PCR. The results are negative.