

HEK293T-M-07e | 305105

HEK293T-M-07e - HEK293T-M-07e

HEK293T-M-07e

Culture Medium	RPMI 1640, w: 2.0 mM β -mercaptoethanol, w: 2.0 g/L NaHCO ₃ (Cytion 820700a)
Supplements	15% FBS, GM-CSF (10 ng/mL), 2.5 g/L HEPES
Doubling time	40 - 46 h
Subculturing	1:6
Fluid renewal	Yes
Freeze medium	10% FBS + 10% DMSO
Thawing and Culturing Cells	1. Thaw cells in a water bath at 37°C. 2. Centrifuge cells at 300 x g for 3 min. 3. Resuspend cells in 10 mL of culture medium. 4. Seed cells into a T25 flask at 70% confluency. 5. Incubate cells at 37°C in 5% CO₂. 6. Harvest cells when they reach 80-90% confluency. 7. Seed cells into a T25 flask at 10% confluency. 8. Incubate cells at 37°C in 5% CO₂.
Incubation Atmosphere	37°C, 5% CO ₂
Flask Coating	Yes

XXXX M-07e | 305105

Freezing Procedure

[illegible]

Shipping Conditions

78°C

Storage Conditions

[illegible]

Sterility

XXXXXXXXXXXXXXXXXXXX XXXX XXXXXXXX XXXXXXXX PCR XXXXXXX XXXXX XXXXXXXXXXXXXXX XXXXXXX XXXXXXXXXXXXXXX

[illegible]