

Cell Line B-LCL-HROC39 (Bc HROC39) | 300863

Cell Line B-LCL-HROC39

Description	B-LCL-HROC39 (Bc HROC39) is a B-cell lymphoma cell line derived from a patient with B-cell lymphoma. It is a continuous cell line that grows in suspension and is characterized by the presence of EBV (Epstein-Barr virus) and the expression of B-cell markers. The cell line is derived from a patient with B-cell lymphoma and is characterized by the presence of EBV (Epstein-Barr virus) and the expression of B-cell markers. The cell line is derived from a patient with B-cell lymphoma and is characterized by the presence of EBV (Epstein-Barr virus) and the expression of B-cell markers. B-LCL-HROC39 (Bc HROC39) is a B-cell lymphoma cell line derived from a patient with B-cell lymphoma. It is a continuous cell line that grows in suspension and is characterized by the presence of EBV (Epstein-Barr virus) and the expression of B-cell markers. The cell line is derived from a patient with B-cell lymphoma and is characterized by the presence of EBV (Epstein-Barr virus) and the expression of B-cell markers. B-LCL-HROC39 (Bc HROC39) is a B-cell lymphoma cell line derived from a patient with B-cell lymphoma. It is a continuous cell line that grows in suspension and is characterized by the presence of EBV (Epstein-Barr virus) and the expression of B-cell markers. The cell line is derived from a patient with B-cell lymphoma and is characterized by the presence of EBV (Epstein-Barr virus) and the expression of B-cell markers.
-------------	--

Organism	Human
Tissue	B-cell lymphoma
Disease	B-cell lymphoma
Synonyms	B-LCL CO39, Bc HROC39

Cell Line B-LCL-HROC39

Age	69 years
Gender	Male
Ethnicity	German
Morphology	Clonal B-cell lymphoma
Cell type	B-cell lymphoma
Growth properties	Continuous cell line

Cell Line B-LCL-HROC39

Citation	B-LCL-HROC39 (Bc HROC39) (Cell Line) Cytion 300863
Biosafety level	2
NCBI_TaxID	9606

CellosaurusAccession CVCL YD56

XXXXXXXXXX XIXX-XXXXXXXXXXXXXXXXXXXX

XXXXXXXXXX

Supplements 10% FBS

Freeze medium α -MEM (Gibco), 10% FBS (Gibco) + 10% DMSO (Gibco) α -MEM (Gibco), 10% FBS (Gibco), 10% DMSO (Gibco), 10% C

Incubation Atmosphere 37°C, 5% CO₂,

XXXX B-LCL-HROC39 (Bc HROC39) | 300863

Flask Coating	Corning® Costar® 384 Well Plates
Freezing Procedure	Seal plates with parafilm and store at -78°C
Shipping Conditions	Ship on dry ice
Storage Conditions	Store at -78°C
Sterility	Use sterile technique