

HEK293T HMY-1 | 305145

HEK293T HMY-1

|                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|-----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description     | HEK293T HMY-1 is a HEK293T cell line expressing the HMY-1 protein. The HMY-1 protein is a member of the HMY family and is involved in the regulation of gene expression. The HMY-1 protein is expressed in HEK293T cells and is secreted into the culture medium. The HMY-1 protein is a dimeric protein and is composed of two subunits. The HMY-1 protein is a member of the HMY family and is involved in the regulation of gene expression. The HMY-1 protein is expressed in HEK293T cells and is secreted into the culture medium. The HMY-1 protein is a dimeric protein and is composed of two subunits. |
| Organism        | HEK293T                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Tissue          | HEK293T                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Disease         | HEK293T                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Metastatic site | HEK293T HEK293T HEK293T HEK293T                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Synonyms        | HMY1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |

HEK293T HMY-1

|                   |                 |
|-------------------|-----------------|
| Age               | 62 HEK293T      |
| Gender            | HEK293T         |
| Morphology        | HEK293T HEK293T |
| Growth properties | HEK293T         |

HEK293T HEK293T HEK293T

|                      |                                       |
|----------------------|---------------------------------------|
| Citation             | HMY-1 (HEK293T HEK293T Cytion 305145) |
| Biosafety level      | 1                                     |
| NCBI_TaxID           | 9606                                  |
| CellosaurusAccession | CVCL_2950                             |

HEK293T HEK293T HEK293T HEK293T

Cytion GmbH | Nikola-Tesla-Str. 3 | 69124 Heidelberg | Heidelberg  
Heidelberg | +49(0)6221 405780 | [www.cytion.com](http://www.cytion.com) | [info@cytion.com](mailto:info@cytion.com)

**HMV-1 | 305145**

**Flask Coating**

Fluoropolymer

**Freezing Procedure**

For freezing, the cells should be seeded into the flask and allowed to attach. The medium should be removed and the cells washed with PBS. The cells should then be trypsinized and resuspended in freezing medium. The suspension should be aliquoted into cryovials and stored at -78°C.

**Shipping Conditions**

The cells should be shipped on dry ice. The shipping container should be pre-cooled and the cells should be placed in the center. The container should be sealed and labeled with the shipping conditions.

**Storage Conditions**

The cells should be stored at -150 to -196 °C in liquid nitrogen. The storage container should be pre-cooled and the cells should be placed in the center. The container should be sealed and labeled with the storage conditions.

**Applications**

**Sterility**

The cells are supplied as a suspension in a sterile medium. The medium is free of mycoplasmas and other contaminants. The cells are suitable for PCR and other applications. The cells should be handled aseptically and the medium should be changed regularly.