

XXX T24 | 300352

XXXXXX XXXXXX

Description	Human T-cell leukemia virus type 1 (HTLV-1) is a retrovirus that causes T-cell hairy-cell leukemia and T-cell hairy-cell lymphoma. It is a member of the Retroviridae family and is classified as a Class II biohazard. HTLV-1 is transmitted through blood, semen, and breast milk. The virus is associated with the development of T-cell hairy-cell leukemia and T-cell hairy-cell lymphoma. HTLV-1 is also associated with the development of T-cell hairy-cell lymphoma. HTLV-1 is a Class II biohazard.
Organism	HTLV-1
Tissue	Human T-cell hairy-cell leukemia and T-cell hairy-cell lymphoma
Disease	T-cell hairy-cell leukemia and T-cell hairy-cell lymphoma
Synonyms	T-24, T 24

XXXXXXXXXXXXXXXXXXXX

Age	82 𐀀𐀀𐀀𐀀
Gender	𐀀𐀀𐀀𐀀
Ethnicity	𐀀𐀀𐀀𐀀𐀀𐀀
Morphology	𐀀𐀀𐀀𐀀 𐀀𐀀𐀀𐀀𐀀
Growth properties	𐀀𐀀𐀀

XIX

Citation	T24 (  Cytion 300352)
Biosafety level	1
NCBI_TaxID	9606
CellosaurusAccession	CVCL_0554

XXXXXXXXXX XXXX-XXXXXXXXXXXXXXX

Antigen expression	HLA A1, A3, B18, Bw35, Cw4, DRw2, Dw4
--------------------	---------------------------------------

Cytion GmbH | Nikola-Tesla-Str. 3 | 69124 Heidelberg | Heidelberg
Heidelberg: +49(0)6221 405780 | www.cytion.com | info@cytion.com

T24 | 300352

**Thawing and
Culturing Cells**

1. 将冻存管放入 37°C 水浴中快速解冻，解冻后立即将细胞悬液转移到含有 1 mL 完全培养基的细胞培养瓶中，轻轻混匀。
2. 将细胞悬液转移到含有 10 mL 完全培养基的细胞培养瓶中，轻轻混匀。将细胞培养瓶放入 37°C 培养箱中培养。
3. 将细胞培养瓶放入 37°C 培养箱中培养，24 小时后，将细胞培养瓶放入 37°C 培养箱中培养。
4. 将细胞培养瓶放入 37°C 培养箱中培养，24 小时后，将细胞培养瓶放入 37°C 培养箱中培养。
5. 将细胞培养瓶放入 37°C 培养箱中培养，24 小时后，将细胞培养瓶放入 37°C 培养箱中培养。
6. 将细胞培养瓶放入 37°C 培养箱中培养，24 小时后，将细胞培养瓶放入 37°C 培养箱中培养。
7. 将细胞培养瓶放入 37°C 培养箱中培养，24 小时后，将细胞培养瓶放入 37°C 培养箱中培养。
8. 将细胞培养瓶放入 37°C 培养箱中培养，24 小时后，将细胞培养瓶放入 37°C 培养箱中培养。

**Incubation
Atmosphere**

37°C, 5% CO₂, 饱和湿度

Flask Coating

无涂层

**Freezing
Procedure**

将细胞悬液转移到含有 1 mL 完全培养基的细胞培养瓶中，轻轻混匀。将细胞培养瓶放入 37°C 培养箱中培养。

**Shipping
Conditions**

将细胞悬液转移到含有 1 mL 完全培养基的细胞培养瓶中，轻轻混匀。将细胞培养瓶放入 37°C 培养箱中培养。

**Storage
Conditions**

将细胞悬液转移到含有 1 mL 完全培养基的细胞培养瓶中，轻轻混匀。将细胞培养瓶放入 37°C 培养箱中培养。

细胞培养

Sterility

将细胞悬液转移到含有 1 mL 完全培养基的细胞培养瓶中，轻轻混匀。将细胞培养瓶放入 37°C 培养箱中培养。

Product sheet

██████T24 | 300352

██████ HLA

A*: 01:01:01

B*: 18:01:01

C*: 05:01:01

DRB1*: 03:01:01

DQA1*: 05:01:01

DQB1*: 02:01:01

DPB1*: 04:01:01

E: 01:01:01