

Product sheet

HEP-CLS-E1 | 400196

HEP-CLS-E1

Description	HEP-CLS-E1, Hepatocellular carcinoma cell line derived from C3H/HE, Hepatocellular carcinoma C3H/HE cell line derived from N-nitrosomethylazoxymethanol (MAM)-induced hepatocellular carcinoma
Organism	HEP
Tissue	HEP
Disease	Hepatocellular carcinoma
Synonyms	HEP-CLS-E1, E1

HEP-CLS-E1

Breed/Subspecies	C3H/HE
Age	HEP
Gender	HEP
Morphology	HEP, HEP
Growth properties	HEP

HEP-CLS-E1

Citation	Hep-CLS-E1 (HEP-CLS-E1 Cytion 400196)
Biosafety level	1
NCBI_TaxID	10090
CellosaurusAccession	CVCL_5777

HEP-CLS-E1

Protein expression	HEP 8, HEP 18, HEP
--------------------	--------------------

Hep-CLS-E1 | 400196

**Thawing and
Culturing Cells**

1. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。
2. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。
3. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。
4. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。
5. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。
6. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。
7. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。
8. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。

**Incubation
Atmosphere**

37°C, 5% CO₂, 饱和湿度

Flask Coating

无涂层

**Freezing
Procedure**

将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。

**Shipping
Conditions**

将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。

**Storage
Conditions**

将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。

细胞培养

Sterility

细胞悬液在37°C CO₂培养箱中培养，用于PCR等下游应用。

细胞悬液在37°C CO₂培养箱中培养，用于PCR等下游应用。