

XXXXXX XXXXXX

XXXXXXXXXXXXXXXXXXXX

XIX

XXXXXXXXXX XIX- XXXXXXXXXXXXXXXX

Product sheet

PK-15 | 607426

Viruses	PCV1 (1) , PCV2 , PCV3
Virus susceptibility	, , , (FMDV), (
Virus resistance	2
Reverse transcriptase	
Culture Medium	EMEM (MEM Eagle), w: 2 mM L-Glutamine, w: 2.2 g/L NaHCO3, w: EBSS (Cytion 820100a)
Supplements	10% FBS 1% NEAA
Dissociation Reagent	
Subculturing	T25, 3-5 ' PBS, 3 3 . , 24 48
Split ratio	1:2 1:4
Seeding density	2×10^4 /
Fluid renewal	2 3
Post-Thaw Recovery	24 48
Freeze medium	(FBS) + 10% DMSO

PK-15 | 607426

Thawing and
Culturing Cells

1. 将细胞悬液缓慢解冻，解冻后将细胞悬液离心，弃去上清液，将细胞重悬于含有生长因子的培养基中。
2. 将细胞悬液接种到含有生长因子的培养基中，在 37°C 下培养 24 小时，使细胞贴壁。
3. 将细胞悬液接种到含有生长因子的培养基中，在 37°C 下培养 24 小时，使细胞贴壁。
4. 将细胞悬液接种到含有生长因子的培养基中，在 37°C 下培养 24 小时，使细胞贴壁。
5. 将细胞悬液接种到含有生长因子的培养基中，在 37°C 下培养 24 小时，使细胞贴壁。
6. 将细胞悬液接种到含有生长因子的培养基中，在 37°C 下培养 24 小时，使细胞贴壁。
7. 将细胞悬液接种到含有生长因子的培养基中，在 37°C 下培养 24 小时，使细胞贴壁。
8. 将细胞悬液接种到含有生长因子的培养基中，在 37°C 下培养 24 小时，使细胞贴壁。

Incubation
Atmosphere

37°C, 5% CO₂, 饱和湿度

Flask Coating

无涂层

Freezing
Procedure

将细胞悬液离心，弃去上清液，将细胞重悬于含有生长因子的培养基中，在 -80°C 下冷冻。

Shipping
Conditions

将细胞悬液离心，弃去上清液，将细胞重悬于含有生长因子的培养基中，在 -80°C 下冷冻。

Storage
Conditions

将细胞悬液离心，弃去上清液，将细胞重悬于含有生长因子的培养基中，在 -150 至 196 K 下冷冻。

Amelogenin STR

Sterility

本产品经过严格的灭菌处理，不含任何细菌、真菌、病毒、支原体、PCR 抑制剂等。本产品不含任何抗生素、血清、生长因子等。

Amelogenin STR

Amelogenin: X,Y