

XXXXXX

Product sheet

HPAF-II | 305088

Incubation
Atmosphere

37°C, 5%_{CO2}, H_2 (5%)

Flask Coating

Yes

Freezing
Procedure

After use, the flask should be stored at -78°C for up to 12 months.

Shipping
Conditions

After use, the flask should be stored at -78°C for up to 12 months.

Storage
Conditions

After use, the flask should be stored at -150 to -196 °C for up to 12 months.

HPAF-II | 305088

Sterility

The flask is sterile and ready for use. PCR products should be stored at -78°C for up to 12 months. The flask is sterile and ready for use. PCR products should be stored at -78°C for up to 12 months.