

[illegible]

333

XXXX XXXXX. XXXXXXX XXX XX

EA. hy 926, EA hy 926, EA-hy926, EAhy 926, EAHY-926, EA.Hy926, EA.hy926, EAhy926, EaHy926, Eahy926

333

XXXXXX

XXX

EA.hy926 (Cytion 305034)

1

9606

CVCL 3901

Cytion GmbH | Nikola-Tesla-Str. 3 | 69124 Heidelberg | Heidelberg
Heidelberg: +49(0)6221 405780 | www.cytion.com | info@cytion.com

EA.hy926 | 305034

**Freezing
Procedure**

Freeze cells in freezing medium (see table) in a 15 ml falcon tube. Freeze at -80°C for 24 hours. Store at -80°C.

**Shipping
Conditions**

Ship cells in freezing medium (see table) in a 15 ml falcon tube. Ship at -80°C.

**Storage
Conditions**

Store cells in freezing medium (see table) in a 15 ml falcon tube at -150 to -196 °C. Store for up to 1 year.

EA.hy926 cells are not certified for clinical use.

Sterility

EA.hy926 cells are not certified for clinical use. PCR genotyping of cells is recommended to confirm identity and purity. Cells are not certified for clinical use. Cells are not certified for clinical use. Cells are not certified for clinical use.