

Product sheet

XXXX XXXX-XXXXXX | 305004

XXXXX XXXXX

Description	Super Tube XXXX XX XXXXX XXXXX XXXXX XX-MDCK, XX XXXXX XX XXX XXXXX XXXXXX XXXX, Canis familiaris. XXXXXXXX XXX XXXXX XXXXXXXX
Organism	XXXX
Tissue	XXXXX XXXXX
Synonyms	MDCK XXXX-XXXXX, SuperTube

XXXXXXXXXXXX

Breed/Subspecies	XXXXX XXXXXXXX
Age	XXXXXX
Gender	XXXXX
Morphology	XXXXXX
Growth properties	XXXX

XXXXXXXXXX XXXXXXXXXXXXXXXX

Citation	Super Tube (XXXXX XXXXXXXX Cytion 305004)
Biosafety level	1
NCBI_TaxID	9615
CellosaurusAccession	CVCL_6453

XXXXXXXXXX XXXX- XXXXXXXXXXXXXXXX

XXXXXXXX

Culture Medium	DMEM:Ham's F12 (1:1), w: 3.1 g/L XXXXXXXX, w: 2.5 mM L-XXXXXXX, w: 15 mM HEPES, w: 0.5 mM XXXX XXXXXXXX, w: 1.2 g/L NaHCO3 820400a)
----------------	---

Cytion GmbH | Nikola-Tesla-Str. 3 | 69124 Heidelberg | Heidelberg
Heidelberg +49(0)6221 405780 | www.cytion.com | info@cytion.com

305004

Freezing
Procedure

Freezing procedure: The sample should be frozen at -78°C.

Shipping
Conditions

Shipping conditions: The sample should be shipped at -78°C.

Storage
Conditions

Storage conditions: The sample should be stored at -150 to -196 °C.

Sterility

The sample is sterile and can be used for PCR. The sample is not suitable for cell culture. The sample is not suitable for sequencing.