

XXXXXX XXXXXX

XXXXXXXXXXXX

XIX XIX XIX XIX XIX XIX XIX XIX XIX XIX

~~XXXXXXXXXX XIX-XXXXXXXXXXXXXXXXXXXX~~

Cytion GmbH | Nikola-Tesla-Str. 3 | 69124 Heidelberg | Heidelberg
Heidelberg +49(0)6221 405780 | www.cytion.com | info@cytion.com

BEAS-2B | 300311

Freezing
Procedure

1. Remove the cells from the culture medium and wash them with PBS. 2. Resuspend the cells in 1 ml of freezing medium. 3. Aliquot the cells into 0.5 ml vials. 4. Freeze the cells at -78°C.

Shipping
Conditions

1. The cells should be shipped on dry ice. 2. The shipping container should be labeled with "Dry Ice" and "Biohazard". 3. The cells should be shipped at -78°C.

Storage
Conditions

1. The cells should be stored at -150 to -196°C. 2. The storage container should be labeled with "Biohazard". 3. The cells should be stored in a liquid nitrogen tank.

BEAS-2B | 300311

Sterility

1. The cells are provided as a frozen sample. 2. The cells are not tested for mycoplasma contamination. 3. The cells are not tested for endotoxin contamination. 4. The cells are not tested for viral contamination. 5. The cells are not tested for bacterial contamination. 6. The cells are not tested for fungal contamination. 7. The cells are not tested for parasitic contamination. 8. The cells are not tested for prion contamination. 9. The cells are not tested for other contaminants.