

SKW-3 | 300343

General information

Description	SKW-3, a human B cell line derived from a 61-year-old male patient with chronic lymphocytic leukemia (CLL), established from peripheral blood mononuclear cells (PBMCs) of patient KE-37, a 61-year-old male patient with acute lymphocytic leukemia (ALL). SKW-3 is a B cell line, while KE-37 is a T cell line. SKW-3 is a B cell line, while KE-37 is a T cell line.
Organism	Human
Tissue	Leukemia
Disease	Chronic lymphocytic leukemia (CLL)
Synonyms	SKW3

Cell characteristics

Age	27 years
Gender	Male
Ethnicity	White
Morphology	Small, round cells
Cell type	B cell
Growth properties	Adherent

Cell line identification

Citation	SKW-3 (SKW3) Cytion 300343
Biosafety level	1
NCBI_TaxID	9606
CellosaurusAccession	CVCL_2197

Product sheet

SKW-3 | 300343

Incubation Atmosphere 37°C, 5% CO₂, humidified

Flask Coating

Freezing Procedure

Shipping Conditions

Storage Conditions

STR

Sterility

- STR
- CSF1PO: 10,12
 - D13S317: 8,12
 - D16S539: 11,12
 - D5S818: 12,13
 - D7S820: 8,12
 - TH01: 6,9,3
 - TPOX: 8
 - vWA: 17,18
 - D3S1358: 15,18
 - D21S11: 28, 29, 39
 - D18S51: 13,18
 - Penta E: 5,14
 - Penta D: 11:15
 - D8S1179: 11,14
 - FGA: 24,25
 - D1S1656: 15.3,16
 - D6S1043: 18,21
 - D2S1338: 19,25
 - D12S391: 19,22
 - D19S433: 13,15

Product sheet

██████████ SKW-3 | 300343

██████████ HLA

A*: '11:01:01, '30:01:01
B*: '35:01:01, '44:02:01
C*: '04:01:01, '05:01:01
DRB1*: '01:03:01, '04:01:01
DQA1*: '01:01:01, '03:03:01
DQB1*: 03:01, 05:01
DPB1*: 04:01:01, 04:02:01
E: 01:01:01