

Cell MKN-45 | 300489

Cell line

Description	Cell MKN-45 is a human gastric cancer cell line, established from a patient with gastric cancer. It is a highly proliferative cell line, suitable for various in vitro studies. Cell MKN-45 is a highly proliferative cell line, suitable for various in vitro studies.
Organism	Human
Tissue	Gastric cancer
Disease	Gastric cancer
Metastatic site	None
Synonyms	MKN 45, MKN45

Cell characteristics

Age	62 years
Gender	Male
Ethnicity	Chinese
Growth properties	Adherent / Suspension

Cell line identification

Citation	MKN-45 (Cell MKN-45) Cytion 300489
Biosafety level	1
NCBI_TaxID	9606
CellosaurusAccession	CVCL_0434

Cell line identification

XXXXXXXXXX

MKN-45 | 300489

**Shipping
Conditions**

Shipping conditions: The product is shipped in a temperature-controlled container at -78°C .

**Storage
Conditions**

Storage conditions: Store at -150 to -196°C in a cryogenic storage container.

Handling

Sterility

The product is sterile and suitable for PCR. It is not recommended to use the product for other applications. The product is not intended for use in clinical or diagnostic applications.