

HOS | 300449

General information

Description	HOS (HOS cell line TE85) is a human osteosarcoma cell line, derived from a 13-year-old female patient with osteosarcoma. HOS cells are highly tumorigenic and are used as a model for studying osteosarcoma. HOS cells are derived from a 143B osteosarcoma cell line.
Organism	Human
Tissue	Cartilage
Disease	Osteosarcoma

Cell characteristics

Age	13 years
Gender	Female
Ethnicity	White
Morphology	Epithelial cells, highly tumorigenic, highly tumorigenic, highly tumorigenic
Growth properties	Highly tumorigenic, highly tumorigenic

Cell line identification

Citation	HOS (HOS cell line Cytion 300449)
Biosafety level	1
NCBI_TaxID	9606
CellosaurusAccession	CVCL_0312

Cell line identification

Isoenzymes	G6PD, B
------------	---------

Cytion GmbH | Nikola-Tesla-Str. 3 | 69124 Heidelberg | Heidelberg
Heidelberg: +49(0)6221 405780 | www.cytion.com | info@cytion.com

37°C, 5% CO₂, [] [] [] [] [] [] [] []

☐ ☐ ☐ ☐ ☐

-78°C

























[illegible]

Abstract

本報刊登廣告，收費低廉，廣告效果顯著，歡迎各界人士垂詢。