

XXXXXX XXXXXX

[illegible]

XXXXXXXXXXXXXXXXXXXX

Growth properties

XIX

CellosaurusAccession CVCL_2580

XXXXXXXXXX XIXX- XXXXXXXXXXXXXXXXXX

Tumorigenic

Products

Karyotype □□□□□□□□, □□□□ □□□□□ = 116, □□□□ □□□□□□□□ □□□□□

XXXXXX

Culture Medium EMEM (MEM Eagle), w: 2 mM L-Glutamine, w: 2.2 g/L NaHCO₃, w: EBSS (XXXX XXXXXX Cytion 820100a)

Supplements □□□□ □□□□□□ 10% FBS □-1% NEAA

Dissociation Reagent 

Subculturing

Seeding density 1 3 x 4 "5

Fluid renewal ☒ ☒ ☒ ☒ ☒ ☒ ☒ ☒

Freeze medium

LMH | 601411

Thawing and
Culturing Cells

1. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。
2. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。
3. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。
4. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。
5. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。
6. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。
7. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。
8. 将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。

Incubation
Atmosphere

37°C, 5% CO₂, 饱和湿度

Flask Coating

无涂层

Freezing
Procedure

将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。

Shipping
Conditions

将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。

Storage
Conditions

将细胞悬液缓慢加入含有培养基的细胞培养瓶中，轻轻混匀，然后将培养瓶放入37°C CO₂培养箱中培养。

细胞培养

Sterility

细胞培养液经过严格的灭菌处理，确保无菌。PCR 试剂和模板 DNA 均经过严格的灭菌处理，确保无菌。