

XXXXXX XXXXXX

XXXXXXXXXX	
Breed/Subspecies	BALB/c
Cell type	XXXXXXXXXX
Growth properties	XXXXX

Citation	WEHI-3B (HEP-2, HEp-2) Cytion 400376)
Biosafety level	2
NCBI_TaxID	10090
CellosaurusAccession	CVCL_2239

Receptors expressed	Viruses
(Fc), (C3)	(HIV receptor) Dengue

Cytion GmbH | Nikola-Tesla-Str. 3 | 69124 Heidelberg | Heidelberg
Heidelberg +49(0)6221 405780 | www.cytion.com | info@cytion.com

WEHI-3B | 400376

Incubation
Atmosphere

37°C, 5%_{CO2}, humidified air

Flask Coating

None

Freezing
Procedure

Cells are harvested by centrifugation at 300g for 5 minutes. The supernatant is removed and the cell pellet is resuspended in freezing medium. The cells are then frozen in a controlled rate freezer and stored at -78°C.

Shipping
Conditions

Cells are shipped in dry ice and stored at -78°C.

Storage
Conditions

Cells are stored in liquid nitrogen at -150 °C for up to 196 days.

Characteristics

Sterility

Cells are tested for mycoplasma contamination by PCR. The results are negative. The cells are also tested for endotoxin contamination and the results are negative.