

**143B Cells | 305232****Informações gerais****Description**

The 143B cell line is a human osteosarcoma cell line derived from a bone tumor. It is often used in cancer research due to its high metastatic potential and ability to form osteosarcoma tumors in vivo. These cells exhibit several key characteristics of osteosarcoma, including the expression of osteoblastic markers and the ability to produce osteoid. The 143B cell line is especially valuable for studying the mechanisms underlying bone cancer progression and metastasis, as well as for testing potential therapeutic agents aimed at treating osteosarcoma.

143B cells are known for their rapid growth and high transfection efficiency, making them suitable for genetic manipulation and various experimental applications. These cells have been utilized in studies investigating the role of specific genes and signaling pathways in osteosarcoma development and resistance to chemotherapy. Additionally, the 143B cell line serves as a model for exploring the interactions between osteosarcoma cells and the bone microenvironment, providing insights into the complex biology of bone tumors.

**Organism** Human**Tissue** Bone, right femur**Disease** Osteosarcoma**Synonyms** 143b, 143 B, 143B TK-, 143B.TK-, 143BTK-, 143TK-, HOS-143B, HOS-143b, GM05887, GM05887A**Características****Age** 13 years**Gender** Female**Ethnicity** Caucasian**Growth properties** Adherent**Dados regulatórios****Citation** 143B (Cytion catalog number 305232)**Biosafety level** 1**NCBI\_TaxID** 9606

**143B Cells | 305232**

CellosaurusAccession CVCL\_2270

**Dados biomoleculares****Manuseio**

**Culture Medium** EMEM (MEM Eagle), w: 2 mM L-Glutamine, w: 2.2 g/L NaHCO<sub>3</sub>, w: EBSS (Cytion article number 820100a)

**Supplements** Supplement the medium with 10% FBS and 1% NEAA

**Dissociation Reagent** Accutase

**Subculturing** Remove the old medium from the adherent cells and wash them with PBS that lacks calcium and magnesium. For T25 flasks, use 3-5 ml of PBS, and for T75 flasks, use 5-10 ml. Then, cover the cells completely with Accutase, using 1-2 ml for T25 flasks and 2.5 ml for T75 flasks. Let the cells incubate at room temperature for 8-10 minutes to detach them. After incubation, gently mix the cells with 10 ml of medium to resuspend them, then centrifuge at 300xg for 3 minutes. Discard the supernatant, resuspend the cells in fresh medium, and transfer them into new flasks that already contain fresh medium.

**Seeding density**  $2 \times 10^4$  cells/cm<sup>2</sup>

**Freeze medium** As a cryopreservation medium, we use complete growth medium (including FBS) + 10% DMSO for adequate post-thaw viability, or CM-1 (Cytion catalog number 800100), which includes optimized osmoprotectants and metabolic stabilizers to enhance recovery and reduce cryo-induced stress.

**143B Cells | 305232**

**Thawing and  
Culturing Cells**

1. Confirm that the vial remains deeply frozen upon delivery, as cells are shipped on dry ice to maintain optimal temperatures during transit.
2. Upon receipt, either store the cryovial immediately at temperatures below -150°C to ensure the preservation of cellular integrity, or proceed to step 3 if immediate culturing is required.
3. For immediate culturing, swiftly thaw the vial by immersing it in a 37°C water bath with clean water and an antimicrobial agent, agitating gently for 40-60 seconds until a small ice clump remains.
4. Perform all subsequent steps under sterile conditions in a flow hood, disinfecting the cryovial with 70% ethanol before opening.
5. Carefully open the disinfected vial and transfer the cell suspension into a 15 ml centrifuge tube containing 8 ml of room-temperature culture medium, mixing gently.
6. Centrifuge the mixture at 300 x g for 3 minutes to separate the cells and carefully discard the supernatant containing residual freezing medium.
7. Gently resuspend the cell pellet in 10 ml of fresh culture medium. For adherent cells, divide the suspension between two T25 culture flasks; for suspension cultures, transfer all the medium into one T25 flask to promote effective cell interaction and growth.
8. Adhere to established subculture protocols for continued growth and maintenance of the cell line, ensuring reliable experimental outcomes.

**Incubation  
Atmosphere**

37°C, 5% CO<sub>2</sub>, humidified atmosphere.

**Shipping  
Conditions**

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately -78 °C throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

**Storage  
Conditions**

For long-term preservation, place vials in vapor-phase liquid nitrogen at about -150 to -196 °C. Storage at -80 °C is acceptable only as a short interim step before transfer to liquid nitrogen.

**Controle de Qualidade e Análise Molecular**

**143B Cells | 305232**

**Sterility**

Mycoplasma contamination is excluded using both PCR-based assays and luminescence-based mycoplasma detection methods.

To ensure there is no bacterial, fungal, or yeast contamination, cell cultures are subjected to daily visual inspections.