

**Wilms1 Cells | 300411**

**Informações gerais**

**Description**

The Wilms1 cell line was derived from a primary Wilms tumor sample obtained from a patient presenting with large bilateral kidney tumors, indicative of Wilms tumor, a pediatric nephroblastoma. This cell line harbors a homozygous nonsense mutation in the WT1 gene (c.149 C>A, p.S50X), which results in a truncated and non-functional WT1 protein. The WT1 gene, critical for kidney development and function, is frequently mutated in Wilms tumor, particularly in those with a stromal subtype that exhibits ectopic mesenchymal differentiation. Wilms1 cells, therefore, represent a unique in vitro model for studying the consequences of WT1 loss of function in tumor biology.

The Wilms1 cell line maintains a stable karyotype with no significant chromosomal abnormalities, allowing for reliable long-term culture. These cells exhibit a mesenchymal phenotype, characterized by the expression of vimentin and the absence of epithelial markers such as cytokeratin, consistent with their stromal origin. Additionally, the cell line demonstrates limited but notable mesenchymal differentiation capacity, including the ability to differentiate into muscle-like cells under appropriate conditions. This makes Wilms1 an invaluable tool for investigating the molecular mechanisms of mesenchymal differentiation and its deregulation in Wilms tumor pathogenesis.

Wilms1 has also been used to study the activation status of key signaling pathways involved in tumor progression. Proteomic analyses have shown that Wilms1 cells exhibit phosphorylation and activation of several receptor tyrosine kinases, including EGFR and PDGFR $\beta$ , as well as downstream MAPK signaling pathways. These findings highlight the relevance of the Wilms1 cell line in exploring targeted therapeutic approaches for Wilms tumor by dissecting the role of these pathways in cancer cell survival, proliferation, and differentiation.

**Organism** Human

**Tissue** Kidney

**Applications** In vitro cell culture model. Biochemical studies

**Synonyms** Wilms1-2l

**Características**

**Age** 2 years

**Gender** Female

**Ethnicity** Caucasian

**Morphology** Spindle-shaped

**Cell type** Wilms cells

**Wilms1 Cells | 300411**

**Growth properties** Adherent

**Dados regulatórios**

**Citation** Wilms1 (Cytion catalog number 300411)

**Biosafety level** 1

**NCBI\_TaxID** 9606

**CellosaurusAccession** CVCL\_A5SC

**Dados biomoleculares**

**Receptors expressed** Receptor tyrosine kinases EGFR, EphA7, PDGFRalpha, FGFR1, PDGFRbeta, Axl

**Tumorigenic** Yes, in nude mice. Forms tumor with small cells consistent with Wilms' tumor (xenografts may not represent Wilm's tumors completely, see E. Kuncz Stroup 2017)

**Viruses** HIV-1: negative, HBV: negative, HCV: negative

**Mutational profile** WT1 mutation status: homozygous c. 149 C>A, p.S50x, LOH: 11p11-11pter, CTNNB1 mutation status: heterozygous TCT>TTT, p.S45F

**Karyotype** 46, normal

**Manuseio**

**Culture Medium** MSCGM kit (from Lonza)

**Dissociation Reagent** Accutase

**Doubling time** 24 hours

## Wilms1 Cells | 300411

**Subculturing** Remove the old medium from the adherent cells and wash them with PBS that lacks calcium and magnesium. For T25 flasks, use 3-5 ml of PBS, and for T75 flasks, use 5-10 ml. Then, cover the cells completely with Accutase, using 1-2 ml for T25 flasks and 2.5 ml for T75 flasks. Let the cells incubate at room temperature for 8-10 minutes to detach them. After incubation, gently mix the cells with 10 ml of medium to resuspend them, then centrifuge at 300xg for 3 minutes. Discard the supernatant, resuspend the cells in fresh medium, and transfer them into new flasks that already contain fresh medium.

**Seeding density**  $1 \times 10^4$  cells/cm<sup>2</sup>

**Fluid renewal** 1 to 2 times per week

**Post-Thaw Recovery** Fast

**Freeze medium** As a cryopreservation medium, we use complete growth medium (including FBS) + 10% DMSO for adequate post-thaw viability, or CM-1 (Cytion catalog number 800100), which includes optimized osmoprotectants and metabolic stabilizers to enhance recovery and reduce cryo-induced stress.

## Wilms1 Cells | 300411

### Thawing and Culturing Cells

1. Confirm that the vial remains deeply frozen upon delivery, as cells are shipped on dry ice to maintain optimal temperatures during transit.
2. Upon receipt, either store the cryovial immediately at temperatures below  $-150^{\circ}\text{C}$  to ensure the preservation of cellular integrity, or proceed to step 3 if immediate culturing is required.
3. For immediate culturing, swiftly thaw the vial by immersing it in a  $37^{\circ}\text{C}$  water bath with clean water and an antimicrobial agent, agitating gently for 40-60 seconds until a small ice clump remains.
4. Perform all subsequent steps under sterile conditions in a flow hood, disinfecting the cryovial with 70% ethanol before opening.
5. Carefully open the disinfected vial and transfer the cell suspension into a 15 ml centrifuge tube containing 8 ml of room-temperature culture medium, mixing gently.
6. Centrifuge the mixture at  $300 \times g$  for 3 minutes to separate the cells and carefully discard the supernatant containing residual freezing medium.
7. Gently resuspend the cell pellet in 10 ml of fresh culture medium. For adherent cells, divide the suspension between two T25 culture flasks; for suspension cultures, transfer all the medium into one T25 flask to promote effective cell interaction and growth.
8. Adhere to established subculture protocols for continued growth and maintenance of the cell line, ensuring reliable experimental outcomes.

### Incubation Atmosphere

$37^{\circ}\text{C}$ , 5%  $\text{CO}_2$ , humidified atmosphere.

### Shipping Conditions

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately  $-78^{\circ}\text{C}$  throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

### Storage Conditions

For long-term preservation, place vials in vapor-phase liquid nitrogen at about  $-150$  to  $-196^{\circ}\text{C}$ . Storage at  $-80^{\circ}\text{C}$  is acceptable only as a short interim step before transfer to liquid nitrogen.

## Controle de Qualidade e Análise Molecular

**Wilms1 Cells | 300411**

**Sterility**

Mycoplasma contamination is excluded using both PCR-based assays and luminescence-based mycoplasma detection methods.

To ensure there is no bacterial, fungal, or yeast contamination, cell cultures are subjected to daily visual inspections.