

**MHH-ES1 Cells | 300136**

**Informações gerais**

**Description**

The MHH-ES1 cell line is derived from a patient with Ewing sarcoma, a highly aggressive bone and soft tissue cancer predominantly affecting children and young adults. This cell line is a valuable model for studying the molecular mechanisms underlying Ewing sarcoma, particularly the role of the EWSR1-FLI1 fusion gene, which is characteristic of this cancer type. The fusion gene results from a translocation between chromosomes 11 and 22, leading to the production of an oncogenic transcription factor that drives tumorigenesis. MHH-ES1, like other Ewing sarcoma cell lines, is utilized to investigate the pathways influenced by EWSR1-FLI1, including alterations in cell proliferation, differentiation, and apoptosis.

Researchers use the MHH-ES1 cell line to evaluate the efficacy of various therapeutic agents targeting pathways critical for Ewing sarcoma survival and proliferation. For example, it is instrumental in testing small molecule inhibitors, RNA interference, and CRISPR-Cas9 gene editing techniques aimed at disrupting the EWSR1-FLI1 fusion gene or its downstream effectors. Additionally, MHH-ES1 serves as a model to study resistance mechanisms to conventional chemotherapy and to identify novel biomarkers for early diagnosis and treatment response monitoring in Ewing sarcoma patients.

**Organism** Human

**Tissue** Bone

**Disease** Ewing's Sarcoma

**Metastatic site** Ascites

**Synonyms** MHH-ES-1, MHHES1

**Características**

**Age** 12 years

**Gender** Male

**Ethnicity** Turkish

**Morphology** Small round cells

**Growth properties** Adherent, clusters

**Dados regulatórios**

## MHH-ES1 Cells | 300136

<b>Citation</b>	MHH-ES1 (Cytion catalog number 300136)
-----------------	--

<b>Biosafety level</b>	1
------------------------	---

<b>NCBI_TaxID</b>	9606
-------------------	------

<b>CellosaurusAccession</b>	CVCL_1411
-----------------------------	-----------

## Dados biomoleculares

## Manuseio

<b>Culture Medium</b>	RPMI 1640, w: 2.0 mM stable Glutamine, w: 2.0 g/L NaHCO <sub>3</sub> (Cytion article number 820700a)
-----------------------	--

<b>Supplements</b>	Supplement the medium with 10% FBS
--------------------	------------------------------------

<b>Dissociation Reagent</b>	Accutase
-----------------------------	----------

<b>Subculturing</b>	Remove the old medium from the adherent cells and wash them with PBS that lacks calcium and magnesium. For T25 flasks, use 3-5 ml of PBS, and for T75 flasks, use 5-10 ml. Then, cover the cells completely with Accutase, using 1-2 ml for T25 flasks and 2.5 ml for T75 flasks. Let the cells incubate at room temperature for 8-10 minutes to detach them. After incubation, gently mix the cells with 10 ml of medium to resuspend them, then centrifuge at 300xg for 3 minutes. Discard the supernatant, resuspend the cells in fresh medium, and transfer them into new flasks that already contain fresh medium.
---------------------	---

<b>Seeding density</b>	1 to 2 x 10 <sup>4</sup> cells/cm <sup>2</sup>
------------------------	--

<b>Fluid renewal</b>	Every 3 to 5 days
----------------------	-------------------

<b>Post-Thaw Recovery</b>	After thawing, plate the cells at 5 x 10 <sup>4</sup> cells/cm <sup>2</sup> and allow the cells to recover from the freezing process and to adhere for at least 24 hours.
---------------------------	---

<b>Freeze medium</b>	As a cryopreservation medium, we use complete growth medium (including FBS) + 10% DMSO for adequate post-thaw viability, or CM-1 (Cytion catalog number 800100), which includes optimized osmoprotectants and metabolic stabilizers to enhance recovery and reduce cryo-induced stress.
----------------------	---

## MHH-ES1 Cells | 300136

### Thawing and Culturing Cells

1. Confirm that the vial remains deeply frozen upon delivery, as cells are shipped on dry ice to maintain optimal temperatures during transit.
2. Upon receipt, either store the cryovial immediately at temperatures below  $-150^{\circ}\text{C}$  to ensure the preservation of cellular integrity, or proceed to step 3 if immediate culturing is required.
3. For immediate culturing, swiftly thaw the vial by immersing it in a  $37^{\circ}\text{C}$  water bath with clean water and an antimicrobial agent, agitating gently for 40-60 seconds until a small ice clump remains.
4. Perform all subsequent steps under sterile conditions in a flow hood, disinfecting the cryovial with 70% ethanol before opening.
5. Carefully open the disinfected vial and transfer the cell suspension into a 15 ml centrifuge tube containing 8 ml of room-temperature culture medium, mixing gently.
6. Centrifuge the mixture at  $300 \times g$  for 3 minutes to separate the cells and carefully discard the supernatant containing residual freezing medium.
7. Gently resuspend the cell pellet in 10 ml of fresh culture medium. For adherent cells, divide the suspension between two T25 culture flasks; for suspension cultures, transfer all the medium into one T25 flask to promote effective cell interaction and growth.
8. Adhere to established subculture protocols for continued growth and maintenance of the cell line, ensuring reliable experimental outcomes.

### Incubation Atmosphere

$37^{\circ}\text{C}$ , 5%  $\text{CO}_2$ , humidified atmosphere.

### Shipping Conditions

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately  $-78^{\circ}\text{C}$  throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

### Storage Conditions

For long-term preservation, place vials in vapor-phase liquid nitrogen at about  $-150$  to  $-196^{\circ}\text{C}$ . Storage at  $-80^{\circ}\text{C}$  is acceptable only as a short interim step before transfer to liquid nitrogen.

## Controle de Qualidade e Análise Molecular

**MHH-ES1 Cells | 300136**

**Sterility**

Mycoplasma contamination is excluded using both PCR-based assays and luminescence-based mycoplasma detection methods.

To ensure there is no bacterial, fungal, or yeast contamination, cell cultures are subjected to daily visual inspections.