

Product sheet

SRA01-04 | 305436

Description	SRA01/04		HLE	
	SRA01/04	130	3.5	
	SRA01/04		---	20%

Organism

Tissue

Synonyms SRA 01-04 HLEC-SRA 01-04 HLE-SRA-01-04 SRA01-04 (HLE) SRA 01/04 HLEC-SRA 01/04 HLE-SRA-01/04 SRA01/04 (HLE)

Age 3

Gender

Morphology

Cell type

Growth properties

Citation SRA01-04 Cytion 305436

Biosafety level 1

NCBI_TaxID 9606

CellosaurusAccession CVCL_7157

GMO Status GMO-S1 SRA01-04 SV40 T

Product sheet

SRA01-04 | 305436

Protein expression	UniProtKB P03070 SV40 T							
Isoenzymes	LD NP							
Culture Medium	DMEM	4.5 g/L	4 mM L-	3.7 g/L NaHCO ₃	1.0 mM	Cytion	820300a	
Supplements	20%		(FBS)					
Dissociation Reagent	Accutase							
Seeding density	2-4*10^4/cm²							
Fluid renewal	2							
Freeze medium	FBS + 10% DMSO				CM-1	Cytion	800100	

Product sheet

SRA01-04 | 305436

Thawing and
Culturing Cells

- | | | | | | |
|----|---------|--------|------|-------|-----|
| 1. | | | | | |
| 2. | | -150°C | | | 3 |
| 3. | | | 37°C | 40-60 | |
| 4. | | | | 70% | |
| 5. | | | 8 | 15 | |
| 6. | 300 x g | 3 | | | |
| 7. | 10 | | | T25 | T25 |
| 8. | | | | | |

Incubation
Atmosphere

37°C 5% CO₂

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 °C -80 °C

Sterility

PCR