

Product sheet

HSF SV40 | 305338

Description	HSF(SV40)	40	SV40	T	T-Ag	SV40 T-Ag	p53
	SV40 T-Ag		HSF(SV40)				

Organism

Metastatic site	SV40
-----------------	------

Applications SV40 p53/Rb

Morphology

Cell type

Growth properties

Citation	HSF(SV40) Cytion 305338
----------	-------------------------

Biosafety level 1

NCBI_TaxID	9606
------------	------

CellosaurusAccession

GMO Status	GMO-S1	HSF	SV40 T
------------	--------	-----	--------

Culture Medium	DMEM	4.5 g/L	4 mM L-	3.7 g/L NaHCO ₃	1.0 mM	Cytion	820300a
----------------	------	---------	---------	----------------------------	--------	--------	---------

Product sheet

HSF SV40 | 305338

Supplements	10%	FBS			
Dissociation Reagent	Accutase				
Doubling time	24	36			
Subculturing	PBS	Accutase	8–10	300×g	3
Split ratio	1	3			
Seeding density	1	3×10^4	/cm ²		
Fluid renewal	2	3			
Freeze medium		FBS + 10% DMSO		CM-1 Cytion	800100
Thawing and Culturing Cells	1.				
	2.	-150°C		3	
	3.	37°C		40-60	
	4.			70%	
	5.	8		15	
	6.	300 x g	3		
	7.	10		T25	T25
	8.				
Incubation Atmosphere	37°C	5% CO ₂			

Product sheet

HSF SV40 | 305338

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 °C

-80 °C

Sterility

PCR