

Product sheet

Hs 611.T | 305915

Description

Hs 611.T 48 cHL T cHL Hs 611.T
Hs 611.T HRS NF-κB JAK-STAT CD30 CD15 CD25)

Organism

Tissue

Disease

Metastatic site

Applications NF-κB JAK-STAT CD30/CD15/CD25

Synonyms Hs-611-T Hs611T HS611T

Age 48

Gender

Ethnicity

Morphology

Cell type B

Growth properties

Citation Hs 611.T Cytion 305915

Biosafety level 1

NCBI_TaxID 9606

Hs 611.T | 305915

CellosaurusAccession CVCL_0823

Culture Medium DMEM 4.5 g/L 4 mM L- 3.7 g/L NaHCO₃ 1.0 mM Cytion 820300a

Supplements 10% FBS

Dissociation Reagent

Doubling time 24 48

Subculturing 3–5 × 10⁵ / 300×g 5

Split ratio 1 3

Seeding density 2 4 × 10⁵

Fluid renewal 2 3

Freeze medium 10% DMSO

Product sheet

Hs 611.T | 305915

Thawing and
Culturing Cells

- | | | | | |
|----|---------|--------|-------|---|
| 1. | | | | |
| 2. | | -150°C | | 3 |
| 3. | | 37°C | 40-60 | |
| 4. | | | 70% | |
| 5. | | 8 | 15 | |
| 6. | 200 x g | 5 | | |
| 7. | | | | |

Incubation
Atmosphere

37°C 5% CO₂

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 °C

-80 °C