

P388-D1 | 400308

Description [P388 D1(IL-1)] -1 IL-1

Organism

Tissue

Disease

Synonyms P-388D1 P388D1 P388.D1 P3 88 D1

Breed/Subspecies DBA/2

Gender

Morphology

Cell type

Growth properties

Citation P388-D1 Cytion 400308

Biosafety level 1

NCBI_TaxID 10090

CellosaurusAccession CVCL_0477

Antigen expression H-2d

P388-D1 | 400308

Tumorigenic

Viruses MAP Ektromelie K 3 PVM LCM MVM GD VII H-1 MHV LDV

Reverse transcriptase

MSI-status

Culture Medium RPMI 1640 2.0 mM 2.0 g/L NaHCO₃ Cytion 820700a

Supplements 10% FBS

Subculturing ⁵ / 5 × 10 3 × 10⁵ 1 × 10⁶ /

Seeding density 1 × 10⁶ /

Fluid renewal

Post-Thaw Recovery ⁶ 24 10

Freeze medium FBS + 10% DMSO CM-1 Cytion 800100

P388-D1 | 400308

**Thawing and
Culturing Cells**

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

**Incubation
Atmosphere**

37°C 5% CO₂

**Shipping
Conditions**

-78 °C

**Storage
Conditions**

-150 -196 °C

-80 °C

Sterility

PCR