

**HK-CRISPR-CAP-H-mEGFP | 301568**

<b>Description</b>	HK-CRISPR-CAP-H-mEGFP	CRISPR-Cas9
	HK-CRISPR-CAP-H-mEGFP	

**Organism**

**Tissue**

**Disease**

**Metastatic site**

<b>Applications</b>	Condensin I	CAP-H	CRISPR
---------------------	-------------	-------	--------

<b>Synonyms</b>	HK-CRISPR-CAP-H-mEGFP #86	HK CRISPR CAP-H-mEGFP
-----------------	---------------------------	-----------------------

<b>Age</b>	30
------------	----

**Gender**

**Ethnicity**

**Morphology**

**Cell type**

**Growth properties**

<b>Citation</b>	HK-CRISPR-CAP-H-mEGFP Cytion 301568
-----------------	-------------------------------------

<b>Biosafety level</b>	1
------------------------	---

<b>NCBI_TaxID</b>	9606
-------------------	------

**HK-CRISPR-CAP-H-mEGFP | 301568**

**CellosaurusAccession** CVCL\_UR43

**Depositor** EMBL

**GMO Status** GMO-S1 HeLa Kyoto CAP-H CRISPR mEGFP

**Products** EGFP

**Culture Medium** DMEM 4.5 g/L 4 mM L- 3.7 g/L NaHCO<sub>3</sub> 1.0 mM Cytion 820300a

**Supplements** 10% FBS

**Dissociation Reagent** Accutase

**Subculturing** PBS T25 3-5 PBS T75 5-10 Accutase T25 1-2 T75

**Fluid renewal** 2 3

**Freeze medium** FBS + 10% DMSO CM-1 Cytion 800100

**HK-CRISPR-CAP-H-mEGFP | 301568**

Thawing and Culturing Cells				
1.				
2.		-150°C		3
3.			37°C	40-60
4.				70%
5.			8	15
6.	300 x g	3		
7.	10		T25	T25
8.				

**Incubation Atmosphere**

37°C 5% CO<sub>2</sub>

Shipping Conditions	
	-78 °C

**Storage Conditions**

-150 -196 °C

-80 °C

Sterility	
	PCR