

NCI-H1650 | 305059

Culture Medium	RPMI 1640	2.0 mM	2.0 g/L NaHCO ₃	Cytion	820700a
-----------------------	-----------	--------	----------------------------	--------	---------

Supplements	10%	FBS
--------------------	-----	-----

Dissociation Reagent	Accutase
-----------------------------	----------

Subculturing	PBS	T25	3-5	PBS	T75	5-10	Accutase	T25	1-2	T75
---------------------	-----	-----	-----	-----	-----	------	----------	-----	-----	-----

Fluid renewal	2	3
----------------------	---	---

Freeze medium	FBS + 10% DMSO	CM-1	Cytion	800100
----------------------	----------------	------	--------	--------

Thawing and Culturing Cells	1.				
	2.	-150°C			3
	3.		37°C		40-60
	4.			70%	
	5.		8	15	
	6.	300 x g	3		
	7.	10		T25	T25
	8.				

Incubation Atmosphere	37°C	5% CO ₂
------------------------------	------	--------------------

NCI-H1650 | 305059

**Shipping
Conditions**

-78 °C

**Storage
Conditions**

-150 -196 °C

-80 °C

Sterility

PCR