

Product sheet

HROC300 T2 M1 | 300866

Description	HROC300 T2 M1	HROC					T2		
	HROC300 T2 M1								
	HROC			CIN	CpG	CIMP	MSS	MSI-H	
Organism									
Tissue									
Disease	TNM	T4aN1bM1R2L0V1	G2	Lk(n) + 3	Σ Lk(n) 22				
Metastatic site	TNM M1	T4aN1bM1R2							
Applications	MSS	HROC Linnebacher							
Age	73								
Gender									
Ethnicity									
Morphology									
Cell type									
Growth properties									
Citation	HROC300 T2 M1	Cytion	300866						
Biosafety level	1								
NCBI_TaxID	9606								

Product sheet

HROC300 T2 M1 | 300866

CellosaurusAccession CVCL_VQ94

GMO Status HROC Linnebacher

MSI-status MSS

Culture Medium DMEM Ham's F12 1:1 3.1 g/L 2.5 mM L- 15 mM HEPES w: 0.5 mM w: 1.2 g/L NaHCO₃ Cytion

Supplements 10% FBS

Dissociation Reagent Accutase

Subculturing PBS T25 3-5 PBS T75 5-10 Accutase T25 1-2 T75

Fluid renewal 3 5

Freeze medium FBS + 10% DMSO CM-1 Cytion 800100

Product sheet

HROC300 T2 M1 | 300866

Thawing and
Culturing Cells

- | | | | | | |
|----|---------|--------|------|-------|-----|
| 1. | | | | | |
| 2. | | -150°C | | | 3 |
| 3. | | | 37°C | 40-60 | |
| 4. | | | | 70% | |
| 5. | | | 8 | 15 | |
| 6. | 300 x g | 3 | | | |
| 7. | 10 | | | T25 | T25 |
| 8. | | | | | |

Incubation
Atmosphere

37°C 5% CO₂

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 °C -80 °C

Sterility

PCR