

EB3 | 300373

**Description**

EB3

EB3

**Organism**

**Tissue**

**Disease**

**Metastatic site**

**Applications**

**Synonyms**

EB-3 - -3 GM04679

**Age**

3

**Gender**

**Ethnicity**

**Morphology**

**Cell type**

B

**Growth properties**

**Citation**

EB3 Cytion 300373

**Biosafety level**

2

**NCBI\_TaxID**

9606

EB3 | 300373

**CellosaurusAccession** CVCL\_1185

**Surface antigens** HLA A3 Aw32 Cw2

**Isoenzymes** G6PD, A

**Viruses** EBV EBNA

**Culture Medium** RPMI 1640 2.0 mM 2.0 g/L NaHCO<sub>3</sub> Cytion 820700a

**Supplements** 10%

**Subculturing** 1 × 10<sup>5</sup> /

**Freeze medium** FBS + 10% DMSO CM-1 Cytion 800100

**Thawing and Culturing Cells**

1.				
2.		-150°C		3
3.		37°C	40-60	
4.			70%	
5.		8	15	
6.	300 x g	3		
7.	10		T25	T25
8.				

