

Product sheet

TERA-C3H | 400360

|             |   |   |     |
|-------------|---|---|-----|
| Description | 1 | 4 | C3H |
|-------------|---|---|-----|

Organism

|        |
|--------|
| Tissue |
|--------|

Disease

|                  |     |
|------------------|-----|
| Breed/Subspecies | C3H |
|------------------|-----|

Age16

|        |
|--------|
| Gender |
|--------|

Morphology

|                   |
|-------------------|
| Growth properties |
|-------------------|

|          |          |        |        |
|----------|----------|--------|--------|
| Citation | TERA-C3H | Cytion | 400360 |
|----------|----------|--------|--------|

Biosafety level1

|            |       |
|------------|-------|
| NCBI_TaxID | 10090 |
|------------|-------|

CellosaurusAccessionCVCL\_AV91

|             |     |
|-------------|-----|
| Tumorigenic | C3H |
|-------------|-----|

VirusesMAPEktromelieKKilhamReo 3PVMLCMMVMDGVIIH-1MHV L

Product sheet

TERA-C3H | 400360

|                             |          |                    |         |                            |         |        |                      |
|-----------------------------|----------|--------------------|---------|----------------------------|---------|--------|----------------------|
| Culture Medium              | DMEM     | 4.5 g/L            | 4 mM L- | 3.7 g/L NaHCO <sub>3</sub> | 1.0 mM  | Cytion | 820300a              |
| Supplements                 | 10%      | FBS                |         |                            |         |        |                      |
| Dissociation Reagent        | Accutase |                    |         |                            |         |        |                      |
| Subculturing                |          | PBS                | T25     | 3-5                        | PBS T75 | 5-10   | Accutase T25 1-2 T75 |
| Fluid renewal               | 2        |                    |         |                            |         |        |                      |
| Freeze medium               |          | FBS + 10% DMSO     |         |                            | CM-1    | Cytion | 800100               |
| Thawing and Culturing Cells | 1.       |                    |         |                            |         |        |                      |
|                             | 2.       |                    | -150°C  |                            |         | 3      |                      |
|                             | 3.       |                    |         | 37°C                       |         | 40-60  |                      |
|                             | 4.       |                    |         |                            | 70%     |        |                      |
|                             | 5.       |                    |         | 8                          | 15      |        |                      |
|                             | 6.       | 300 x g            | 3       |                            |         |        |                      |
|                             | 7.       | 10                 |         |                            | T25     |        | T25                  |
|                             | 8.       |                    |         |                            |         |        |                      |
| Incubation Atmosphere       | 37°C     | 5% CO <sub>2</sub> |         |                            |         |        |                      |
| Shipping Conditions         |          |                    |         |                            |         | -78 °C |                      |
| Storage Conditions          |          |                    | -150    | -196 °C                    |         | -80 °C |                      |

Product sheet

TERA-C3H | 400360

Sterility

PCR