

P3HR1 | 302310

Description P3HR1 - Jiyoye Jiyoye 7

Organism

Disease

Applications B HLA

Synonyms P 3 HR 1 P3HR-1 P3HRI P3HR1-BL P3J HR-1 P3J-HR-1 P3JHR-1 P3JHR1 P3J HR1-K P3J.HR1K HR-1 HR1K PO

Age 7

Gender

Morphology

Cell type B

Growth properties

Citation P3HR1 Cytion 302310

Biosafety level 2

NCBI_TaxID 9606

CellosaurusAccession CVCL_2676

Antigen expression CD10+ CD19+ CD20+ CD21- CD22+ CD23- CD24- CD37+ CD32+ CD38+ CD39- CD40+ CD54+ CD71+ CD72+ CD73- CD74- CD75+ CD77- MHC II -

P3HR1 | 302310

Karyotype 46

Culture Medium RPMI 1640 2.0 mM 2.0 g/L NaHCO₃ Cytion 820700a

Supplements 10% FBS

Subculturing ⁵ / 5 × 10⁵ 3 × 10⁵ 1 × 10⁶ /

Seeding density 3 × 10⁵ /

Post-Thaw Recovery

Freeze medium FBS + 10% DMSO CM-1 Cytion 800100

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

Incubation Atmosphere 37°C 5% CO₂

P3HR1 | 302310

**Shipping
Conditions**

-78 °C

**Storage
Conditions**

-150 -196 °C

-80 °C

Sterility

PCR