

Product sheet

Ba/F3 | 305224

Description	BA/F3	BALB/c	B	BaF3					IL-3
	Baf3	IL-3		Ba/F3	F3	Ba/F3	F3		
		BA/F3				CML	BCR-ABL	Ba/F3	
		BA/F3							

Organism

Tissue	
--------	--

Synonyms BA/F3, BaF3, BAF3, Baf3

Breed/Subspecies	C3H
------------------	-----

Morphology

Cell type	B
-----------	---

Growth properties

Citation	Ba/F3 Cytion 305224
----------	---------------------

Biosafety level 1

NCBI_TaxID	10090
------------	-------

CellosaurusAccession CVCL_0161

Karyotype	Ba/F3	33%
-----------	-------	-----

Ba/F3 | 305224

Culture Medium	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO ₃	Cytion	820700a
-----------------------	--------------------	------------------------------	--------	---------

Supplements	5%	FBS	10 ng/mL	IL-3
--------------------	----	-----	----------	------

Subculturing	⁵ / ⁵ 1×10 ⁶ / 3×10	5×10
---------------------	--	------

Freeze medium	FBS +10% DMSO	CM-1	Cytion	800100	CM-1	Cytion	800100
----------------------	---------------	------	--------	--------	------	--------	--------

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

Incubation Atmosphere	37°C, 5% CO ₂
------------------------------	--------------------------

Flask Coating

Freezing Procedure	-78 °C
---------------------------	--------

Shipping Conditions	-78 °C
----------------------------	--------

Product sheet

Ba/F3 | 305224

Storage Conditions -150 -196 -80 °C

Sterility PCR