

Product sheet

HEK293-CTLA4 | 305423

Description	/ 6,250					
	HEK293-CTLA4	HEK293	CTLA4	4,000	inscreenex	CTLA4
		CTLA4				

Organism

Tissue						
--------	--	--	--	--	--	--

Disease	/	HEK293				
---------	---	--------	--	--	--	--

Applications	CTLA-4	T				
--------------	--------	---	--	--	--	--

Age						
-----	--	--	--	--	--	--

Gender

Morphology						
------------	--	--	--	--	--	--

Cell type

Growth properties						
-------------------	--	--	--	--	--	--

Citation	HEK293-CTLA4	Cytion	305423			
----------	--------------	--------	--------	--	--	--

Biosafety level

1

NCBI_TaxID	9606					
------------	------	--	--	--	--	--

CellosaurusAccession

CVCL_6G20

HEK293-CTLA4 | 305423

GMO Status	GMO-S1	HEK293	CTLA-4
-------------------	--------	--------	--------

Receptors expressed	CTLA4 (CD152)
----------------------------	---------------

Culture Medium	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO ₃	Cytion	820700a
-----------------------	--------------------	------------------------------	--------	---------

Supplements	10% FBS 1 mM	10 mM HEPES	1% NEAA	G418-	1 mg/mL
--------------------	--------------	-------------	---------	-------	---------

Dissociation Reagent	-EDTA
-----------------------------	-------

Doubling time	24-36
----------------------	-------

Subculturing	PBS	PBS5%	2-3	/EDTA	T25	1	T75	3	37°C
---------------------	-----	-------	-----	-------	-----	---	-----	---	------

Split ratio	1 5
--------------------	-----

Seeding density	2 4×10 ⁴ /
------------------------	-----------------------

Fluid renewal	2 3
----------------------	-----

Post-Thaw Recovery	1:2 1:3	T25	24
---------------------------	---------	-----	----

Freeze medium	FBS +10% DMSO	CM-1 Cytion	800100	CM-1 Cytion	800100
----------------------	---------------	-------------	--------	-------------	--------

Product sheet

HEK293-CTLA4 | 305423

Thawing and
Culturing Cells

- | | | | | | |
|----|---------|--------|-----|-------|-----|
| 1. | | | | | |
| 2. | | -150°C | | | 3 |
| 3. | | 37°C | | 40-60 | |
| 4. | | | 70% | | |
| 5. | | 8 | 15 | | |
| 6. | 300 x g | 3 | | | |
| 7. | 10 | | | T25 | T25 |
| 8. | | | | | |

Incubation
Atmosphere

37°C, 5%_{CO2}

Flask Coating

Freezing
Procedure

-78 °C

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 -80 °C

Sterility

PCR