

Product sheet

HEK293-FOLR1 | 305425

Description	/ 6,250				
	HEK293-FOLR1	HEK293	FOLR1	15,000	inscreenex “ ”
		FOLR1			

Organism

Tissue					
--------	--	--	--	--	--

Disease	/	HEK293			
---------	---	--------	--	--	--

Applications	FOLR1	ADC	ADCC/CDC		
--------------	-------	-----	----------	--	--

Age					
-----	--	--	--	--	--

Gender

Morphology					
------------	--	--	--	--	--

Cell type

Growth properties					
-------------------	--	--	--	--	--

Citation	HEK293-FOLR1	Cytion	305425		
----------	--------------	--------	--------	--	--

Biosafety level

1

NCBI_TaxID	9606				
------------	------	--	--	--	--

CellosaurusAccession

CVCL_6G24

HEK293-FOLR1 | 305425

GMO Status	GMO-S1	HEK293	FOLR1
-------------------	--------	--------	-------

Receptors expressed	FOLR1	α FR α	FBP
----------------------------	-------	----------------------	-----

Culture Medium	RPMI 1640	w 2.0 mM	w 2.0 g/L NaHCO ₃	Cytion	820700a
-----------------------	-----------	----------	------------------------------	--------	---------

Supplements	10% FBS	1 mM	10 mM HEPES	1% NEAA	G418-	1 mg/mL
--------------------	---------	------	-------------	---------	-------	---------

Dissociation Reagent	-EDTA
-----------------------------	-------

Doubling time	24-36
----------------------	-------

Subculturing	PBS	PBS5%	2-3	/EDTA	T25	1	T75	3	37°C
---------------------	-----	-------	-----	-------	-----	---	-----	---	------

Split ratio	1	5
--------------------	---	---

Seeding density	2×10^4	/
------------------------	-----------------	---

Fluid renewal	2	3
----------------------	---	---

Post-Thaw Recovery	1:2	1:3	T25	24
---------------------------	-----	-----	-----	----

Freeze medium	FBS +10% DMSO	CM-1	Cytion	800100	CM-1	Cytion	800100
----------------------	---------------	------	--------	--------	------	--------	--------

Product sheet

HEK293-FOLR1 | 305425

Thawing and
Culturing Cells

1.							
2.		-150°C				3	
3.		37°C			40-60		
4.			70%				
5.			8	15			
6.	300 x g	3					
7.	10				T25		T25
8.							

Incubation
Atmosphere

37°C, 5%_{CO2}

Flask Coating

Freezing
Procedure

-78 °C

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 -80 °C

Sterility

PCR