

Product sheet

HSF (SV40) | 305338

Description	HSF SV40	Simian Virus 40 (SV40) T (T-Ag) SV40 T-Ag	p53
	HSF SV40	SV40 T-Ag	

Organism

Metastatic site	SV40
-----------------	------

Applications SV40 p53/Rb

Morphology

Cell type

Growth properties

Citation	HSF SV40 Cytion 305338
----------	------------------------

Biosafety level 1

NCBI_TaxID	9606
------------	------

CellosaurusAccession

GMO Status	GMO-S1 HSF SV40 T
------------	-------------------

Culture Medium	DMEM:Ham's F12 1:1 w: 3.1 g/L w: 2.5 mM L- w: 15 mM HEPES w: 0.5 mM w: 1.2 g/L NaHCO3 Cytion
----------------	--

Product sheet

HSF (SV40) | 305338

Supplements 10% FBS 50 /

Dissociation Reagent Accutase

Doubling time 24 36

Subculturing PBS Accutase 8–10 300×g 3

Split ratio 1 3

Seeding density 1 3 × 10⁴ /cm²

Fluid renewal 2 3

Freeze medium FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

Incubation Atmosphere 37°C, 5% CO₂

Flask Coating

Product sheet

HSF (SV40) | 305338

Freezing
Procedure

-78 °C

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 -80 °C

Sterility

PCR